


## ONCOLOGY INFUSION SET INSTRUCTIONS FOR USE



### RELATED PRODUCTS

One – way, two – way and multi – way oncology infusion sets

### PRODUCT DESCRIPTION

They are sets that allow oncology treatments to be administered with the help of gravity or appropriate infusion pump through a single closed set.

### INSTRUCTIONS FOR USE

1. Close the clamps of infusion set.
2. Open the protective cap of spike tip. Plunge the spike into the port of saline or similar solution bag to be used as a washing solution in one stroke.
3. Open the clamp near the spike tip. Allow the drip chamber fill halfway with liquid and allow the liquid move through the set.
4. Remove the air inside the set by opening all the clamps and by applying pressure to the bag containing the liquid. Close all clamps.
5. Prepare the transfer set in accordance with the instructions for use. (See Oncology Transfer Set User Manual).
6. Connect the distal end of infusion set to the patient's IV device.
7. Lock the luer tip of the properly prepared transfer set into one of the check valves/needle-free connectors of the infusion set.
8. Check the patient's vascular access by opening all the clamps on the infusion set. Close the clamp near the spike tip of the infusion set, leave the other clamps open.
9. Perform the drug infusion by opening the clamp of the transfer set.

10. Close the clamp of the transfer set when the infusion is finished. Open all clamps of the infusion set and perform washing with the solution which is connected.
11. If different drugs are to be administered, repeat the same procedure for different drugs over the check valves/needle-free connectors of the two-way or multiple-way infusion set.
12. Close all clamps of the infusion set when the procedures are completed. Dispose of the set as a closed system without separating its parts.

### WARNINGS & PRECAUTIONS

1. It should be used by healthcare personnel.
2. NOT SUITABLE for blood products containing cells.
3. DO NOT USE if the package has been opened or damaged.
4. DO NOT USE expired products.
5. Check that all protective caps are present before use.
6. DO NOT ACCESS needle-free valves with a needle or blunt cannula.
7. DO NOT USE in patients who are hypersensitive to any of the materials used.
8. If the silicone part of the valve becomes damaged or does not return to its original form after multiple access, replace the set. This may cause leakage.
9. DO NOT OVERTIGHTEN the connections.
10. If it is leaking, replace the product with a new one.
11. Replace the product with a new one in case the set separates from the bonding areas.
12. The product is for single use only. DO NOT RESTERILIZE.
13. DO NOT LEAVE infusion sets unattended.

## ONCOLOGY INFUSION SET INSTRUCTIONS FOR USE

14. If you observe that a drug or liquid has damaged the infusion set, replace the set with a new one and inform our company by naming the drug.
15. When the drug administration process is finished, dispose of the infusion set and the connected transfer sets as a whole without separating them.
16. Follow local regulations regarding the disposal of hazardous drugs.

### STORAGE TEMPERATURE

Product storage temperature is 15 - 30°C.

### TIME OF USE

Product usage time is 48 hours.



















### SHELF LIFE

The shelf life of the product is 5 years.

### USERS

Specialist physicians, nurses, pharmacists and other relevant health-trained workers.

### LABEL DISCLOSURES

				
Manufacturer Information	Notified Body Number	Manufacturing Date	Expiry Date	Lot Number
				
Reference Number	Sterilized with EO Gas	Do Not Re-sterilize	Do Not Use If Package is Damaged	Keep Away From The Sun
				
See User Manual	Storage Temperature Limits	Single Use Only	Caution	Pyrogen Free
				
DEHP – Free	Latex – Free	Medical Device		