





EVA

BAGS



EVA BAG



EVA

BAGS

EVA INFUSION BAGS w. FLL/MLL LUER CONNECTOR – 03

EVA INFUSION BAGS w. BIG BORE CONNECTOR – 04

EVA INFUSION BAGS w. NEEDLE - FREE CONNECTOR – 05

EVA INFUSION BAGS w. EXTENSION LINE – 06

EVA INFUSION BAGS w. 3 - WAY TRANSFER SET – 07



EVA INFUSION BAGS w. 6 - WAY TRANSFER SET – 09



EVA INFUSION BAGS w. FLL/MLL LUER CONNECTOR

Ethylene vinyl acetate bags with male or female luer ends. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options are available.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.



TRANSPARENT EVA BAGS w. FLL LUER CONNECTOR	
REF NO	VOLUME (mL)
OE250	250
OE500	500
OE1000	1000
OE2000	2000
OE3000	3000
OE5000	5000

UV PROTECTED EVA BAGS w. FLL LUER CONNECTOR	
REF NO	VOLUME (mL)
OE250P	250
OE500P	500
OE1000P	1000
OE2000P	2000
OE3000P	3000
OE5000P	5000

TRANSPARENT EVA BAGS w. MLL LUER CONNECTOR	
REF NO	VOLUME (mL)
OE251	250
OE501	500
OE1001	1000
OE2001	2000
OE3001	3000
OE5001	5000



UV PROTECTED EVA BAGS w. MLL LUER CONNECTOR	
REF NO	VOLUME (mL)
OE251P	250
OE501P	500
OE1001P	1000
OE2001P	2000
OE3001P	3000
OE5001P	5000



EVA INFUSION BAGS w. BIG BORE CONNECTOR

Ethylene vinyl acetate bags with male or female luer ends, containing big bore connectors. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.



TRANSPARENT EVA BAGS w. BIG BORE & FLL CONNECTOR	
REF NO	VOLUME (mL)
OE250B	250
OE500B	500
OE1000B	1000
OE2000B	2000
OE3000B	3000
OE5000B	5000

UV PROTECTED EVA BAGS w. BIG BORE & FLL CONNECTOR	
REF NO	VOLUME (mL)
OE250PB	250
OE500PB	500
OE1000PB	1000
OE2000PB	2000
OE3000PB	3000
OE5000PB	5000

TRANSPARENT EVA BAGS w. BIG BORE & MILL CONNECTOR	
REF NO	VOLUME (mL)
OE251B	250
OE501B	500
OE1001B	1000
OE2001B	2000
OE3001B	3000
OE5001B	5000


UV PROTECTED EVA BAGS w. BIG BORE & MILL CONNECTOR	
REF NO	VOLUME (mL)
OE251PB	250
OE501PB	500
OE1001PB	1000
OE2001PB	2000
OE3001PB	3000
OE5001PB	5000



EVA INFUSION BAGS w. NEEDLE - FREE CONNECTOR

Ethylene vinyl acetate bags containing a needle-free connector at the port end. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.



TRANSPARENT EVA BAGS w. NEEDLE-FREE CONNECTOR	
REF NO	VOLUME (mL)
OE250V	250
OE500V	500
OE1000V	1000
OE2000V	2000
OE3000V	3000
OE5000V	5000


UV PROTECTED EVA BAGS w. NEEDLE-FREE CONNECTOR	
REF NO	VOLUME (mL)
OE250PV	250
OE500PV	500
OE1000PV	1000
OE2000PV	2000
OE3000PV	3000
OE5000PV	5000



EVA INFUSION BAGS w. EXTENSION LINE

Ethylene vinyl acetate bags containing an extension line. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.



TRANSPARENT EVA BAGS w. EXTENSION LINE & FLL CONNECTOR	
REF NO	VOLUME (mL)
OE250-11	250
OE500-11	500
OE1000-11	1000
OE2000-11	2000
OE3000-11	3000
OE5000-11	5000

TRANSPARENT EVA BAGS w. EXTENSION LINE & BIG BORE CONNECTOR	
REF NO	VOLUME (mL)
OE250B-11	250
OE500B-11	500
OE1000B-11	1000
OE2000B-11	2000
OE3000B-11	3000
OE5000B-11	5000

UV PROTECTED EVA BAGS w. EXTENSION LINE & FLL CONNECTOR	
REF NO	VOLUME (mL)
OE250P-12	250
OE500P-12	500
OE1000P-12	1000
OE2000P-12	2000
OE3000P-12	3000
OE5000P-12	5000


UV PROTECTED EVA BAGS w. EXTENSION LINE & BIG BORE CONNECTOR	
REF NO	VOLUME (mL)
OE250PB-12	250
OE500PB-12	500
OE1000PB-12	1000
OE2000PB-12	2000
OE3000PB-12	3000
OE5000PB-12	5000



EVA INFUSION BAGS w. 3 - WAY TRANSFER SET

Ethylene vinyl acetate bags designed for separate storage of lipid emulsion or other special solutions. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.





EVA INFUSION BAGS w. 3 - WAY TRANSFER SET

TRANSPARENT EVA BAGS w. 3-WAY TRANSFER SET & FLL CONNECTOR

REF NO	VOLUME (mL)
OE250-30	250
OE500-30	500
OE1000-30	1000
OE2000-30	2000
OE3000-30	3000
OE5000-30	5000

UV PROTECTED EVA BAGS w. 3-WAY TRANSFER SET & FLL CONNECTOR

REF NO	VOLUME (mL)
OE250P-30	250
OE500P-30	500
OE1000P-30	1000
OE2000P-30	2000
OE3000P-30	3000
OE5000P-30	5000

TRANSPARENT EVA BAGS w. 3-WAY TRANSFER SET & MLL CONNECTOR

REF NO	VOLUME (mL)
OE251-30	250
OE501-30	500
OE1001-30	1000
OE2001-30	2000
OE3001-30	3000
OE5001-30	5000

UV PROTECTED EVA BAGS w. 3-WAY TRANSFER SET & MLL CONNECTOR

REF NO	VOLUME (mL)
OE251P-30	250
OE501P-30	500
OE1001P-30	1000
OE2001P-30	2000
OE3001P-30	3000
OE5001P-30	5000

TRANSPARENT EVA BAGS w. 3-WAY TRANSFER SET & BIG BORE CONNECTOR

REF NO	VOLUME (mL)
OE250B-30	250
OE500B-30	500
OE1000B-30	1000
OE2000B-30	2000
OE3000B-30	3000
OE5000B-30	5000

UV PROTECTED EVA BAGS w. 3-WAY TRANSFER SET & BIG BORE CONNECTOR

REF NO	VOLUME (mL)
OE250PB-30	250
OE500PB-30	500
OE1000PB-30	1000
OE2000PB-30	2000
OE3000PB-30	3000
OE5000PB-30	5000

TRANSPARENT EVA BAGS w. 3-WAY TRANSFER SET & NEEDLE-FREE CONNECTOR

REF NO	VOLUME (mL)
OE250V-30	250
OE500V-30	500
OE1000V-30	1000
OE2000V-30	2000
OE3000V-30	3000
OE5000V-30	5000

UV PROTECTED EVA BAGS w. 3-WAY TRANSFER SET & NEEDLE-FREE CONNECTOR

REF NO	VOLUME (mL)
OE250PV-30	250
OE500PV-30	500
OE1000PV-30	1000
OE2000PV-30	2000
OE3000PV-30	3000
OE5000PV-30	5000



EVA INFUSION BAGS w. 6 - WAY TRANSFER SET

Ethylene vinyl acetate bags designed for separated storage of lipid emulsion or other special solutions. UV protected options are available. Available in different volumes from 1000mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 1.000 mL - 2.000 mL - 3.000 mL
5.000 mL volume options.



TRANSPARENT EVA BAGS w. 6-WAY TRANSFER SET & FLL CONNECTOR

REF NO	VOLUME (mL)
OE1000-40	1000
OE2000-40	2000
OE3000-40	3000
OE5000-40	5000

TRANSPARENT EVA BAGS w. 6-WAY TRANSFER SET & MLL CONNECTOR

REF NO	VOLUME (mL)
OE1001-40	1000
OE2001-40	2000
OE3001-40	3000
OE5001-40	5000

TRANSPARENT EVA BAGS w. 6-WAY TRANSFER SET & BIG BORE CONNECTOR

REF NO	VOLUME (mL)
OE1000B-40	1000
OE2000B-40	2000
OE3000B-40	3000
OE5000B-40	5000

IUV PROTECTED EVA BAGS w. 6-WAY TRANSFER SET & FLL CONNECTOR

REF NO	VOLUME (mL)
OE1000P-40	1000
OE2000P-40	2000
OE3000P-40	3000
OE5000P-40	5000

IUV PROTECTED EVA BAGS w. 6-WAY TRANSFER SET & MLL CONNECTOR

REF NO	VOLUME (mL)
OE1001P-40	1000
OE2001P-40	2000
OE3001P-40	3000
OE5001P-40	5000

IUV PROTECTED EVA BAGS w. 6-WAY TRANSFER SET & BIG BORE CONNECTOR

REF NO	VOLUME (mL)
OE1000PB-40	1000
OE2000PB-40	2000
OE3000PB-40	3000
OE5000PB-40	5000

ONCOERA
Safe for your health

ERASER[®]
MEDICAL

ERASER MEDICAL LLC

KEMALPAŞA MH. PINAR CD. NO:78
PINARBAŞI/BORNOVA, İZMİR - TURKEY

TEL : +90 232 479 80 60

FAX : +90 232 479 80 25

www.eraser.com.tr
info@eraser.com.tr
info@oncoera.com.tr

