




**NUTRITION**

**SETS**



## NUTRITION

## SETS

**ENTERAL NUTRITION SETS – 03**

**ENTERAL ADMINISTRATION SETS – 05**

**ENTERAL EXTENSION SETS – 06**

**NEONATAL PARENTERAL NUTRITION SETS – 07**

PEDIATRIC SETS w. 0,2 MICRON FILTER – **08**

PEDIATRIC SETS w. 1,2 MICRON FILTER – **09**

PEDIATRIC SETS w. 1,2 MICRON FILTER - 160 CM – **10**


TWO - WAY PEDIATRIC SETS w. 0,2 MICRON FILTER – **11**

TWO - WAY PEDIATRIC SETS w. 1,2 MICRON FILTER – **12**

**ADULT PARENTERAL NUTRITION SETS – 13**

ADULT SETS w. 0,2 MICRON FILTER – **14**


ADULT SETS w. 1,2 MICRON FILTER – **14**



## ENTERAL NUTRITION SETS

NUTRIERA Enteral Nutrition Sets are designed for feeding patients who cannot take nutrients orally. It has a transparent bag structure where the nutritional liquid can be seen easily, an administration line for transferring the liquid, and a label area where patient information can be written.

- DEHP and LATEX free.
- Includes Enfit connector end with a cap.
- Includes luer lock connector compatible with PEG connectors/jejunal probes.
- 500mL - 800mL - 1000mL - 1200mL - 1500mL - 1800mL - 2000mL bag options.



## ENTERAL NUTRITION SETS



### ENTERAL NUTRITION SETS

SETS USED w. FEEDING PUMP	VOLUME (mL)
EN210P1 Enteral Nutrition Set	500
EN220P1 Enteral Nutrition Set	800
EN230P1 Enteral Nutrition Set	1000
EN240P1 Enteral Nutrition Set	1200
EN250P1 Enteral Nutrition Set	1500
EN260P1 Enteral Nutrition Set	1800
EN270P1 Enteral Nutrition Set	2000

### ENTERAL NUTRITION SETS - DOUBLE BAG


SETS USED w. FEEDING PUMP	VOLUME (mL)
EN234P1 Enteral Nutrition Set w. Flush Bag	1000
EN244P1 Enteral Nutrition Set w. Flush Bag	1200
EN254P1 Enteral Nutrition Set w. Flush Bag	1500

### ENTERAL NUTRITION SETS

SETS USED w. GRAVITY	VOLUME (mL)
EN210G1 Enteral Nutrition Set	500
EN220G1 Enteral Nutrition Set	800
EN230G1 Enteral Nutrition Set	1000
EN240G1 Enteral Nutrition Set	1200
EN250G1 Enteral Nutrition Set	1500
EN260G1 Enteral Nutrition Set	1800
EN270G1 Enteral Nutrition Set	2000

### ENTERAL NUTRITION SETS - BOTTLE CAPPED

ENTERAL NUTRITION SETS - BOTTLE CAPPED
EN185P1 Enteral Nutrition Set w. Bottle Cap



## ENTERAL ADMINISTRATION SETS

NUTRIERA Enteral Administration Sets are designed for feeding patients who cannot take nutrients orally. Sets are used by connecting to enteral bags. Contains an administration line for liquid transfer and a label where the patient information will be written.

- DEHP and LATEX free.
- Includes Enfit connector end with a cap.
- Includes luer lock connector compatible with PEG connectors/jejunal probes.
- Available to use with enteral feeding pump or gravity flow.

### SETS USED w. FEEDING PUMP

EN180P1 Enteral Administration Set

### SETS USED w. GRAVITY


EN180G1 Enteral Administration Set




## ENTERAL EXTENSION SETS

NUTRIERA Enteral Extension Sets are available in varying lengths according to user needs. It can be connected to the Y connector of enteral nutrition sets and enteral administration sets.

- DEHP and LATEX free.
- Includes Enfit connector end with a cap.
- Includes luer lock connector compatible with PEG connectors/jejunal probes.




ENTERAL EXTENSION SETS	LENGTH (cm)
EN301U1 Enteral Extension Set	15
EN401U1 Enteral Extension Set	30
EN501U1 Enteral Extension Set	60
EN601U1 Enteral Extension Set	90
EN701U1 Enteral Extension Set	100
EN801U1 Enteral Extension Set	120




## NEONATAL PARENTERAL NUTRITION SETS

NUTRIERA Neonatal Parenteral Nutrition Sets are configured in accordance with the principles of TPN administration in pediatrics. Contains 0,2  $\mu\text{m}$  or 1,2  $\mu\text{m}$  filter. Allows for proper filtering of nutritional content and lipid content separately to support two-way administration.

- Provides removal of air bubbles in liquid infusion.
- Removes large particles, bacteria and endotoxins from the system.
- Prevents harmful effects of daylight with its UV light-protected structure.
- Supports patient safety.



## NEONATAL PARENTERAL NUTRITION SETS PEDIATRIC SETS w. 0,2 MICRON FILTER




### PVC SETS (DEHP - FREE)

REF NO	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV100P1.102	✓	-	1,5X2,5	✓
IV100S1.102	✓	-	1,5X2,5	-
IV101P1.102	✓	-	1,5X2,7	✓
IV101S1.102	✓	-	1,5X2,7	-
IV102P1.102	✓	-	1,5X3,0	✓
IV102S1.102	✓	-	1,5X3,0	-
IV110P1.102	✓	✓	1,5X2,5	✓
IV110S1.102	✓	✓	1,5X2,5	-
IV111P1.102	✓	✓	1,5X2,7	✓
IV111S1.102	✓	✓	1,5X2,7	-
IV112P1.102	✓	✓	1,5X3,0	✓
IV112S1.102	✓	✓	1,5X3,0	-


### POLYURETHANE SETS

REF NO	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV100P0.102	✓	-	1,5X2,5	✓
IV100S0.102	✓	-	1,5X2,5	-
IV101P0.102	✓	-	1,5X2,7	✓
IV101S0.102	✓	-	1,5X2,7	-
IV102P0.102	✓	-	1,5X3,0	✓
IV102S0.102	✓	-	1,5X3,0	-
IV110P0.102	✓	✓	1,5X2,5	✓
IV110S0.102	✓	✓	1,5X2,5	-
IV111P0.102	✓	✓	1,5X2,7	✓
IV111S0.102	✓	✓	1,5X2,7	-
IV112P0.102	✓	✓	1,5X3,0	✓
IV112S0.102	✓	✓	1,5X3,0	-





## NEONATAL PARENTERAL NUTRITION SETS PEDIATRIC SETS w. 1,2 MICRON FILTER




### PVC SETS (DEHP - FREE)

REF NO	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV100P1.112	✓	-	1,5X2,5	✓
IV100S1.112	✓	-	1,5X2,5	-
IV101P1.112	✓	-	1,5X2,7	✓
IV101S1.112	✓	-	1,5X2,7	-
IV102P1.112	✓	-	1,5X3,0	✓
IV102S1.112	✓	-	1,5X3,0	-
IV110P1.112	✓	✓	1,5X2,5	✓
IV110S1.112	✓	✓	1,5X2,5	-
IV111P1.112	✓	✓	1,5X2,7	✓
IV111S1.112	✓	✓	1,5X2,7	-
IV112P1.112	✓	✓	1,5X3,0	✓
IV112S1.112	✓	✓	1,5X3,0	-

### POLYURETHANE SETS

REF NO	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV100P0.112	✓	-	1,5X2,5	✓
IV100S0.112	✓	-	1,5X2,5	-
IV101P0.112	✓	-	1,5X2,7	✓
IV101S0.112	✓	-	1,5X2,7	-
IV102P0.112	✓	-	1,5X3,0	✓
IV102S0.112	✓	-	1,5X3,0	-
IV110P0.112	✓	✓	0,65X2,5	✓
IV110S0.112	✓	✓	0,65X2,5	-
IV111P0.112	✓	✓	1,5X2,7	✓
IV111S0.112	✓	✓	1,5X2,7	-
IV112P0.112	✓	✓	1,5X3,0	✓
IV112S0.112	✓	✓	1,5X3,0	-

## NEONATAL PARENTERAL NUTRITION SETS PEDIATRIC SETS w. 1,2 MICRON FILTER - 160 CM




### PVC SETS (DEHP - FREE)

REF NO	1,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV100P1.112-160	✓	1,5X2,5	✓
IV100S1.112-160	✓	1,5X2,5	-
IV101P1.112-160	✓	1,5X2,7	✓
IV101S1.112-160	✓	1,5X2,7	-
IV102P1.112-160	✓	1,5X3,0	✓
IV102S1.112-160	✓	1,5X3,0	-
IV104P1.112-160	✓	0,65X2,5	✓
IV104S1.112-160	✓	0,65X2,5	-


### POLYURETHANE SETS

REF NO	1,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV100P0.112-160	✓	1,5X2,5	✓
IV100S0.112-160	✓	1,5X2,5	-
IV101P0.112-160	✓	1,5X2,7	✓
IV101S0.112-160	✓	1,5X2,7	-
IV102P0.112-160	✓	1,5X3,0	✓
IV102S0.112-160	✓	1,5X3,0	-
IV104P0.112-160	✓	0,65X2,5	✓
IV104S0.112-160	✓	0,65X2,5	-



## NEONATAL PARENTERAL NUTRITION SETS

### TWO - WAY PEDIATRIC SETS w. 0,2 MICRON FILTER




#### PVC SETS (DEHP - FREE)

REF NO	0,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV400P1.102	✓	1,5X2,5	✓
IV400S1.102	✓	1,5X2,5	-
IV401P1.102	✓	1,5X2,7	✓
IV401S1.102	✓	1,5X2,7	-
IV402P1.102	✓	1,5X3,0	✓
IV402S1.102	✓	1,5X3,0	-

#### POLYURETHANE SETS

REF NO	0,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV400P0.102	✓	1,5X2,5	✓
IV400S0.102	✓	1,5X2,5	-
IV401P0.102	✓	1,5X2,7	✓
IV401S0.102	✓	1,5X2,7	-
IV402P0.102	✓	1,5X3,0	✓
IV402S0.102	✓	1,5X3,0	-

## NEONATAL PARENTERAL NUTRITION SETS TWO - WAY PEDIATRIC SETS w. 1,2 MICRON FILTER




### PVC SETS (DEHP - FREE)

REF NO	1,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV400P1.112	✓	1,5X2,5	✓
IV400S1.112	✓	1,5X2,5	-
IV401P1.112	✓	1,5X2,7	✓
IV401S1.112	✓	1,5X2,7	-
IV402P1.112	✓	1,5X3,0	✓
IV402S1.112	✓	1,5X3,0	-

### POLYURETHANE SETS


REF NO	1,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV400P0.112	✓	1,5X2,5	✓
IV400S0.112	✓	1,5X2,5	-
IV401P0.112	✓	1,5X2,7	✓
IV401S0.112	✓	1,5X2,7	-
IV402P0.112	✓	1,5X3,0	✓
IV402S0.112	✓	1,5X3,0	-



## ADULT PARENTERAL NUTRITION SETS

NUTRIERA Adult Parenteral Nutrition Sets are configured in accordance with the principles of TPN administration in adults. Contains 0,2  $\mu\text{m}$  or 1,2  $\mu\text{m}$  filter. Allows for proper filtering of nutritional content and lipid content separately to support two-way administration.

- Provides removal of air bubbles in liquid infusion.
- Removes large particles, bacteria and endotoxins from the system.
- Prevents the harmful effects of daylight with its UV light-protected structure.
- Supports patient safety.



## ADULT PARENTERAL NUTRITION SETS ADULT SETS w. 0,2 MICRON FILTER



### PVC SETS (DEHP - FREE)

REF NO	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV103P1.002	✓	-	3,0X4,1	✓
IV103S1.002	✓	-	3,0X4,1	-
IV113P1.002	✓	-	3,0X4,1	✓
IV113S1.002	✓	-	3,0X4,1	-

### POLYURETHANE SETS

REF NO	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV103P0.002	✓	-	3,0X4,1	✓
IV103S0.002	✓	-	3,0X4,1	-
IV113P0.002	✓	-	3,0X4,1	✓
IV113S0.002	✓	-	3,0X4,1	-

## ADULT PARENTERAL NUTRITION SETS ADULT SETS w. 1,2 MICRON FILTER



### PVC SETS (DEHP - FREE)

REF NO	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV103P1.012	✓	-	3,0X4,1	✓
IV103S1.012	✓	-	3,0X4,1	-
IV113P1.012	✓	-	3,0X4,1	✓
IV113S1.012	✓	-	3,0X4,1	-

### POLYURETHANE SETS

REF NO	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV103P0.012	✓	-	3,0X4,1	✓
IV103S0.012	✓	-	3,0X4,1	-
IV113P0.012	✓	-	3,0X4,1	✓
IV113S0.012	✓	-	3,0X4,1	-



**ONCOERA**  
*Safe for your health*

**ERASER**<sup>®</sup>  
MEDICAL

**ERASER MEDICAL LLC**

KEMALPAŞA MH. PINAR CD. NO:78  
PINARBAŞI/BORNOVA, İZMİR - TURKEY

**TEL :** +90 232 479 80 60

**FAX :** +90 232 479 80 25

[www.eraser.com.tr](http://www.eraser.com.tr)  
[info@eraser.com.tr](mailto:info@eraser.com.tr)  
[info@oncoera.com.tr](mailto:info@oncoera.com.tr)

