

ONCOERA

Safe for your health

ONCOERA

Safe for your health

CHAPTERS

COMPANY
PROFILE

PRODUCTION
FACILITY

ONCOLOGY

EVA
BAG

NUTRITION

INFUSION

OEM

COMPOUNDING
DEVICES

PRODUCT CATALOGUE

ERASER[®]
MEDICAL

ORICO

INDEX

COMPANY PROFILE
01

PRODUCTION FACILITY
02

ONCOLOGY
07



EVA BAG
27

NUTRITION
36

INFUSION
50

OEM
79

COMPOUNDING DEVICES
81



ALL IN SAFE WITH
ERASER[®]





ONCOERA



Eraser Medical LLC was founded in 2007 in İzmir. The company registered its experience with the ONCOERA brand, which it has been producing since 2014. Today, Eraser Medical is the manufacturer, authorized dealer and exporter of ONCOERA, NUTRIERA, ONCOMIX and ERAVERSE branded products.

Eraser Medical LLC, which produces disposable sterile & non-sterile medical supplies, also manufactures aseptic drug preparation devices and software, by its strong R&D infrastructure.

Company has adopted a scientific-based and customer-oriented production and sales policy. Products of Eraser Medical are exported to 4 continents.



PRODUCTION FACILITY



The production facility has a closed area of 1,700 m². The production, assembly and packaging of medical devices is carried out in 2 clean rooms of Class 10,000 and 2 clean rooms of Class 100,000. The biological load on the products is kept under control in these clean rooms and sterilization safety is ensured. Disposable medical equipment production consists of injection, extrusion, assembly, packaging and sterilization sections. The storage processes of the products before loading are monitored in a 500 m² area. The raw materials used in production are biocompatible and meet USP and ISO 10933 standards. Eraser Medical LLC has ISO 9001, ISO 13485 and CE certificates.

In the production facility, ONCOERA, NUTRIERA, ONCOMIX and ERAVERSE branded products registered to Eraser Medical LLC are produced.

The R&D center is structured to develop software and device hardware and has a wide range of device park. R&D processes are advanced by channeling the years of experience of the marketing, engineering, software, medical and quality teams to a common point in order to develop technology. In research and development activities, both the needs of the sector are taken into account and the most up-to-date scientific trends in the field of health are utilized.





PRODUCTION FACILITY







PRODUCTION FACILITY



ERASER[®]
MEDİKAL

ONCOERA
Safe for your health

nutriera

ERAVERVE
MEDICAL DEVICES



ONCOLOGY

PRODUCTS



ONCOLOGY

PRODUCTS

VIAL ADAPTERS – 09

SYRINGE ADAPTERS – 10

ONCOLOGY TRANSFER SETS – 11

NEEDLE - FREE CONNECTORS (VALVES) – 12

ONCOLOGY INFUSION SETS – 13

ONE - WAY ONCOLOGY INFUSION SETS – 14

TWO - WAY ONCOLOGY INFUSION SETS – 16

MULTI - WAY ONCOLOGY INFUSION SETS – 18

PACLITAXEL INFUSION SETS – 20

ONCOLOGY INFUSION ADAPTERS – 22

ONCOLOGY PRODUCTS USER MANUAL – 23



VIAL ADAPTERS



ONCOERA Vial Adapters are equipped with a 0,2 µm hydrophobic air filter to protect users from aerosols and drugs from microbiological contamination.

- Minimizes toxic aerosols in the air.
- Neutralizes vial pressure.
- Increases safety by minimizing the risk of surface contamination.
- Protects compounded drug, personnel, and compounding environment.
- Maintains the sterile integrity of drugs.
- Reduces the risk of infection and cross-contamination.
- Provides low dead volume and high flow rate.



VA101F.13



VA101F.20



VA102F.U



VA112F.20



VA104F.20



VA105F.U



VA111F.20



SYRINGE ADAPTERS



ONCOERA SAFE Needle-Free Syringe Adapters are connectors used for the safe preparation and transport of hazardous drugs. Syringe adapters creates a closed system between syringe and vial adapters, as well as administration sets and needle-free connectors.

- Supports closed system.
- Protects compounded drug, personnel, and compounding environment.
- Reduces surface contamination and prevents leakage.
- Quickly and easily connects to a needle-free vein valve.
- Reduces the risk of personnel exposure to hazardous drugs.
- Prevents needle stick accidents.
- Provides safe preparation, handling and administration.
- Compatible with all standard luer connectors.
- Complies with NIOSH requirements for transporting of hazardous drugs.
- Chemical compatibility has been proven.
- Suitable for use with IV administration sets.
- Provides high flow rate.



SF201



SF401



ONCOLOGY TRANSFER SETS



ONCOERA LINE Oncology Transfer Sets are extension lines with needle-free connectors, which allow the connection of prepared drugs to the administration sets and ensure the continuity of the closed system. Transfer sets may include an in-line filter or a needle-free connector to meet different requirements. It can be designed with or without UV protection and set lengths can be configured as desired.

- Supports closed system.
- Protects compounded drug, personnel, and compounding environment.
- UV protected set options are available for drugs that require UV protection.
- Set options with 0,22 µm filters are available for drugs such as paclitaxel.
- Integrated valve options are available to prevent free flow.
- DEHP-Free.

PVC SETS (DEHP - FREE)


REF NO	NEEDLE - FREE CONNECTOR	0.2 MICRON FILTER	UV PROTECTION	LENGTH
LNT001P1	-	-	✓	25 CM
LNT001S1	-	-	-	25 CM
LNT003P1	-	✓	✓	25 CM
LNT003S1	-	✓	-	25 CM
LNT101P1	✓	-	✓	22 CM
LNT101S1	✓	-	-	22 CM
LNT103P1	✓	✓	✓	22 CM
LNT103S1	✓	✓	-	22 CM

POLYURETHANE SETS

REF NO	NEEDLE - FREE CONNECTOR	0.2 MICRON FILTER	UV PROTECTION	LENGTH
LNT001P2	-	-	✓	25 CM
LNT001S2	-	-	-	25 CM
LNT003P2	-	✓	✓	25 CM
LNT003S2	-	✓	-	25 CM
LNT101P2	✓	-	✓	22 CM
LNT101S2	✓	-	-	22 CM
LNT103P2	✓	✓	✓	22 CM
LNT103S2	✓	✓	-	22 CM



NEEDLE - FREE CONNECTORS (VALVES)



ONCOERA IV Needle-Free Connectors are the part of the split septum valve group. Needle-free connectors have become an integral part of intravenous therapy in recent years. The most important reason for this widespread acceptance is reduces the risk of infection of the patient and healthcare personnel. The CDC also recommends the use of split septum vein valves in intravenous therapy to reduce the risk of infection. Needle-free connectors create a closed system with syringe adapters to improve employee and environmental safety.

- Supports closed system.
- Protects compounded drug, personnel, and compounding environment.
- Suitable to use with all luer-lock connectors.
- Forms a microbiological barrier for 7 days.
- Easy to disinfect.
- Prevents needle stick accidents.
- Provides low dead volume and high flow rate.
- Compatible with blood products and oncology drugs.
- DEHP-Free.



IV201
VALVE w. TRANSPARENT SILICONE




IV301
VALVE w. BLUE SILICONE



ONCOLOGY INFUSION SETS


Antineoplastic drugs used in cancer treatment should be prepared, administered and disposed under special conditions due to their toxic effects. Administration of compounded drugs at the right speed and at the right dose is extremely important for the effectiveness of treatment. Oncology infusion sets are the safest solution used in chemotherapy administration and prioritizing patient comfort. Sets contain needle-free connectors that ensure the continuity of closed system.

- Supports closed system.
- Protects compounded drug, personnel, and environment.
- UV protected set options are available.
- Set options with 0,22 µm filters are available, for drugs such as paclitaxel.
- Integrated valve options are available to prevent free flow.
- Removes particles from the system with a 15 µm particle filter in drip chamber.
- Multiple infusions of fluids can be done quickly and safely through a single set.
- Reduces workload significantly and saves time and cost.
- DEHP-Free.






ONCOLOGY INFUSION SETS ONE - WAY ONCOLOGY INFUSION SETS





ONCOLOGY INFUSION SETS ONE - WAY ONCOLOGY INFUSION SETS



PVC SETS (DEHP - FREE)

REF NO	SPIKE	CLAMP	CHECK VALVE	DRIP CHAMBER	ROLLER CLAMP	SILICONE LINE	STOPCOCK	IN - LINE CHECK VALVE	NEEDLE - FREE CONNECTOR	MALE LUER END	CHECK VALVE END	UV PROTECTION	LENGTH
LN113P1	✓	✓	1 PC	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	✓	195CM
LN113S1	✓	✓	1 PC	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	-	195CM
LN151P1	✓	✓	1 PC	✓	✓	24CM SILICONE	-	-	✓	✓	-	✓	220CM
LN151S1	✓	✓	1 PC	✓	✓	24CM SILICONE	-	-	✓	✓	-	-	220CM
LN341P1	✓	✓	1 PC	✓	✓	-	-	-	✓	-	✓	✓	220CM
LN341S1	✓	✓	1 PC	✓	✓	-	-	-	✓	-	✓	-	220CM
LN351P1	✓	✓	1 PC	✓	✓	-	-	-	✓	✓	-	✓	220CM
LN351S1	✓	✓	1 PC	✓	✓	-	-	-	✓	✓	-	-	220CM

Silicone line optimizes the sensitivity of the administration set to deliver fluid through the infusion pump. Gravity infusion sets do not have a silicone line.

POLYURETHANE SETS

REF NO	SPIKE	CLAMP	CHECK VALVE	DRIP CHAMBER	ROLLER CLAMP	SILICONE LINE	STOPCOCK	IN - LINE CHECK VALVE	NEEDLE - FREE CONNECTOR	MALE LUER END	CHECK VALVE END	UV PROTECTION	LENGTH
LN113P2	✓	✓	1 PC	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	✓	195CM
LN113S2	✓	✓	1 PC	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	-	195CM
LN151P2	✓	✓	1 PC	✓	✓	24CM SILICONE	-	-	✓	✓	-	✓	220CM
LN151S2	✓	✓	1 PC	✓	✓	24CM SILICONE	-	-	✓	✓	-	-	220CM
LN341P2	✓	✓	1 PC	✓	✓	-	-	-	✓	-	✓	✓	220CM
LN341S2	✓	✓	1 PC	✓	✓	-	-	-	✓	-	✓	-	220CM
LN351P2	✓	✓	1 PC	✓	✓	-	-	-	✓	✓	-	✓	220CM
LN351S2	✓	✓	1 PC	✓	✓	-	-	-	✓	✓	-	-	220CM

Silicone line optimizes the sensitivity of the administration set to deliver fluid through the infusion pump. Gravity infusion sets do not have a silicone line.




ONCOLOGY INFUSION SETS TWO - WAY ONCOLOGY INFUSION SETS





ONCOLOGY INFUSION SETS TWO - WAY ONCOLOGY INFUSION SETS



PVC SETS (DEHP - FREE)

REF NO	SPIKE	CLAMP	CHECK VALVE	DRIP CHAMBER	ROLLER CLAMP	SILICONE LINE	STOPOCK	IN - LINE CHECK VALVE	NEEDLE - FREE CONNECTOR	MALE LUER END	CHECK VALVE END	UV PROTECTION	LENGTH
LN114P1	✓	✓	2 PCS	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	✓	195CM
LN114S1	✓	✓	2 PCS	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	-	195CM
LN152P1	✓	✓	2 PCS	✓	✓	24CM SILICONE	-	-	✓	✓	-	✓	220CM
LN152S1	✓	✓	2 PCS	✓	✓	24CM SILICONE	-	-	✓	✓	-	-	220CM
LN342P1	✓	✓	2 PCS	✓	✓	-	-	-	✓	-	✓	✓	220CM
LN342S1	✓	✓	2 PCS	✓	✓	-	-	-	✓	-	✓	-	220CM
LN352P1	✓	✓	2 PCS	✓	✓	-	-	-	✓	✓	-	✓	220CM
LN352S1	✓	✓	2 PCS	✓	✓	-	-	-	✓	✓	-	-	220CM

Silicone line optimizes the sensitivity of the administration set to deliver fluid through the infusion pump. Gravity infusion sets do not have a silicone line.


POLYURETHANE SETS

REF NO	SPIKE	CLAMP	CHECK VALVE	DRIP CHAMBER	ROLLER CLAMP	SILICONE LINE	STOPOCK	IN - LINE CHECK VALVE	NEEDLE - FREE CONNECTOR	MALE LUER END	CHECK VALVE END	UV PROTECTION	LENGTH
LN114P2	✓	✓	2 PCS	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	✓	195CM
LN114S2	✓	✓	2 PCS	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	-	195CM
LN152P2	✓	✓	2 PCS	✓	✓	24CM SILICONE	-	-	✓	✓	-	✓	220CM
LN152S2	✓	✓	2 PCS	✓	✓	24CM SILICONE	-	-	✓	✓	-	-	220CM
LN342P2	✓	✓	2 PCS	✓	✓	-	-	-	✓	-	✓	✓	220CM
LN342S2	✓	✓	2 PCS	✓	✓	-	-	-	✓	-	✓	-	220CM
LN352P2	✓	✓	2 PCS	✓	✓	-	-	-	✓	✓	-	✓	220CM
LN352S2	✓	✓	2 PCS	✓	✓	-	-	-	✓	✓	-	-	220CM

Silicone line optimizes the sensitivity of the administration set to deliver fluid through the infusion pump. Gravity infusion sets do not have a silicone line.




ONCOLOGY INFUSION SETS MULTI - WAY ONCOLOGY INFUSION SETS





ONCOLOGY INFUSION SETS MULTI - WAY ONCOLOGY INFUSION SETS



PVC SETS (DEHP - FREE)

REF NO	SPIKE	CLAMP	CHECK VALVE	DRIP CHAMBER	ROLLER CLAMP	SILICONE LINE	STOPCOCK	IN - LINE CHECK VALVE	NEEDLE - FREE CONNECTOR	MALE LUER END	CHECK VALVE END	UV PROTECTION	LENGTH
LN110P1	✓	✓	4 PCS	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	✓	198CM
LN110S1	✓	✓	4 PCS	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	-	198CM
LN154P1	✓	✓	4 PCS	✓	✓	24CM SILICONE	-	-	✓	✓	-	✓	223CM
LN154S1	✓	✓	4 PCS	✓	✓	24CM SILICONE	-	-	✓	✓	-	-	223CM
LN344P1	✓	✓	4 PCS	✓	✓	-	-	-	✓	-	✓	✓	240CM
LN344S1	✓	✓	4 PCS	✓	✓	-	-	-	✓	-	✓	-	240CM
LN354P1	✓	✓	4 PCS	✓	✓	-	-	-	✓	✓	-	✓	240CM
LN354S1	✓	✓	4 PCS	✓	✓	-	-	-	✓	✓	-	-	240CM

Silicone line optimizes the sensitivity of the administration set to deliver fluid through the infusion pump. Gravity infusion sets do not have a silicone line.


POLYURETHANE SETS

REF NO	SPIKE	CLAMP	CHECK VALVE	DRIP CHAMBER	ROLLER CLAMP	SILICONE LINE	STOPCOCK	IN - LINE CHECK VALVE	NEEDLE - FREE CONNECTOR	MALE LUER END	CHECK VALVE END	UV PROTECTION	LENGTH
LN110P2	✓	✓	4 PCS	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	✓	198CM
LN110S2	✓	✓	4 PCS	✓	✓	24CM SILICONE	✓	✓	✓	✓	-	-	198CM
LN154P2	✓	✓	4 PCS	✓	✓	24CM SILICONE	-	-	✓	✓	-	✓	223CM
LN154S2	✓	✓	4 PCS	✓	✓	24CM SILICONE	-	-	✓	✓	-	-	223CM
LN344P2	✓	✓	4 PCS	✓	✓	-	-	-	✓	-	✓	✓	240CM
LN344S2	✓	✓	4 PCS	✓	✓	-	-	-	✓	-	✓	-	240CM
LN354P2	✓	✓	4 PCS	✓	✓	-	-	-	✓	✓	-	✓	240CM
LN354S2	✓	✓	4 PCS	✓	✓	-	-	-	✓	✓	-	-	240CM

Silicone line optimizes the sensitivity of the administration set to deliver fluid through the infusion pump. Gravity infusion sets do not have a silicone line.




ONCOLOGY INFUSION SETS PACLITAXEL INFUSION SETS





ONCOLOGY INFUSION SETS PACLITAXEL INFUSION SETS



PVC SETS (DEHP - FREE)


REF NO	DRIP CHAMBER w. SPIKE	ROLLER CLAMP	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	Y INJECTION PORT	MALE LUER END	VENTED CAP	UV PROTECTION	LENGTH
LN001P1	✓	✓	✓	-	-	✓	✓	✓	180CM
LN001S1	✓	✓	✓	-	-	✓	✓	-	180CM
LN002P1	✓	✓	✓	-	-	✓	✓	✓	200CM
LN002S1	✓	✓	✓	-	-	✓	✓	-	170CM
LN102P1	✓	✓	✓	✓	-	✓	✓	✓	170CM
LN102S1	✓	✓	✓	✓	-	✓	✓	-	170CM
LN121P1	✓	✓	✓	-	✓	✓	✓	✓	180CM
LN121S1	✓	✓	✓	-	✓	✓	✓	-	180CM

POLYURETHANE SETS

REF NO	DRIP CHAMBER w. SPIKE	ROLLER CLAMP	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	Y INJECTION PORT	MALE LUER END	VENTED CAP	UV PROTECTION	LENGTH
LN001P2	✓	✓	✓	-	-	✓	✓	✓	180CM
LN001S2	✓	✓	✓	-	-	✓	✓	-	180CM
LN002P2	✓	✓	✓	-	-	✓	✓	✓	200CM
LN002S2	✓	✓	✓	-	-	✓	✓	-	170CM
LN102P2	✓	✓	✓	✓	-	✓	✓	✓	170CM
LN102S2	✓	✓	✓	✓	-	✓	✓	-	170CM
LN121P2	✓	✓	✓	-	✓	✓	✓	✓	180CM
LN121S2	✓	✓	✓	-	✓	✓	✓	-	180CM



ONCOLOGY INFUSION SETS ONCOLOGY INFUSION ADAPTERS



PVC ADAPTERS (DEHP - FREE)

REF NO	SPIKE	CLAMP	NEEDLE - FREE CONNECTOR	CHECK VALVE	MEMBRANE CONNECTOR	UV PROTECTION	LENGTH
LN104P1	✓	✓	✓	-	✓	✓	12CM
LN105P1	✓	✓	-	2	✓	✓	24CM
LN106P1	✓	✓	-	4	✓	✓	21CM

POLYURETHANE ADAPTERS

REF NO	SPIKE	CLAMP	NEEDLE - FREE CONNECTOR	CHECK VALVE	MEMBRANE CONNECTOR	UV PROTECTION	LENGTH
LN104P2	✓	✓	✓	-	✓	✓	12CM
LN105P2	✓	✓	-	2	✓	✓	24CM
LN106P2	✓	✓	-	4	✓	✓	21CM



ONCOLOGY PRODUCTS USER MANUAL



- 1
- Close the clamp of transfer set. Attach the transfer set to spike port of the IV bag and fill the set with fluid by opening the clamp. Close the clamp again.



- 2
- Open the sterile package of the syringe adapter and syringe.



- 3
- Lock the syringe adapter onto luer lock syringe by twisting.



- 4
- Remove the vial adapter from its sterile package and open the protective cap of vial.



ONCOLOGY PRODUCTS USER MANUAL

5



- Plunge spike of the vial adapter onto the vial stopper in one stroke.

6



- Take the syringe with syringe adapter attached and the vial with vial adapter attached.

7



- Lock the syringe adapter to vial adapter by twisting to provide a closed system.

8



- Withdraw the required amount of liquid into the syringe.
(During withdrawing, the syringe should be held at a 90-degree angle to the ground.)



ONCOLOGY PRODUCTS USER MANUAL



- Open the cap of needle-free connector on the transfer set.



- Lock the syringe with syringe adapter attached to needle-free connector of the transfer set by twisting.



- Transfer required amount of liquid from the syringe into the bag by pushing the syringe plunger. After fluid transfer is complete, withdraw a small amount of fluid from the IV bag into the syringe and return it to the bag. (Flushing)



- Bring prepared medication bag near the infusion set.



ONCOLOGY PRODUCTS USER MANUAL

13



- Lock the luer end of the transfer set to one of the check valves on the infusion set by twisting. Once this connection is established, it should not be opened again.

14



- Open the clamp of transfer set just before drug administration begins. Close the clamps when administration is finished. Dispose of the entire set as cytotoxic waste without disconnecting the infusion set and transfer set.

IMPORTANT

- Aseptic working principles should be followed at all stages of the hazardous drug preparation and administration process. Personal protective equipment should be used appropriately.
- Some medications require special preparation conditions. Follow the manufacturer's instructions.
- This manual contains general usage information. Please contact for detailed information such as reconstitution.

info@oncoera.com.tr ; info@eraser.com.tr





EVA

BAGS



EVA BAG



EVA

BAGS

EVA INFUSION BAGS w. FLL/MLL LUER CONNECTOR – 29

EVA INFUSION BAGS w. BIG BORE CONNECTOR – 30

EVA INFUSION BAGS w. NEEDLE - FREE CONNECTOR – 31

EVA INFUSION BAGS w. EXTENSION LINE – 32

EVA INFUSION BAGS w. 3 - WAY TRANSFER SET – 33



EVA INFUSION BAGS w. 6 - WAY TRANSFER SET – 35



EVA INFUSION BAGS w. FLL/MLL LUER CONNECTOR

Ethylene vinyl acetate bags with male or female luer ends. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options are available.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.



TRANSPARENT EVA BAGS w. FLL LUER CONNECTOR	
REF NO	VOLUME (mL)
OE250	250
OE500	500
OE1000	1000
OE2000	2000
OE3000	3000
OE5000	5000

UV PROTECTED EVA BAGS w. FLL LUER CONNECTOR	
REF NO	VOLUME (mL)
OE250P	250
OE500P	500
OE1000P	1000
OE2000P	2000
OE3000P	3000
OE5000P	5000

TRANSPARENT EVA BAGS w. MLL LUER CONNECTOR	
REF NO	VOLUME (mL)
OE251	250
OE501	500
OE1001	1000
OE2001	2000
OE3001	3000
OE5001	5000



UV PROTECTED EVA BAGS w. MLL LUER CONNECTOR	
REF NO	VOLUME (mL)
OE251P	250
OE501P	500
OE1001P	1000
OE2001P	2000
OE3001P	3000
OE5001P	5000



EVA INFUSION BAGS w. BIG BORE CONNECTOR

Ethylene vinyl acetate bags with male or female luer ends, containing big bore connectors. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.



TRANSPARENT EVA BAGS w. BIG BORE & FLL CONNECTOR	
REF NO	VOLUME (mL)
OE250B	250
OE500B	500
OE1000B	1000
OE2000B	2000
OE3000B	3000
OE5000B	5000

UV PROTECTED EVA BAGS w. BIG BORE & FLL CONNECTOR	
REF NO	VOLUME (mL)
OE250PB	250
OE500PB	500
OE1000PB	1000
OE2000PB	2000
OE3000PB	3000
OE5000PB	5000

TRANSPARENT EVA BAGS w. BIG BORE & MILL CONNECTOR	
REF NO	VOLUME (mL)
OE251B	250
OE501B	500
OE1001B	1000
OE2001B	2000
OE3001B	3000
OE5001B	5000


UV PROTECTED EVA BAGS w. BIG BORE & MILL CONNECTOR	
REF NO	VOLUME (mL)
OE251PB	250
OE501PB	500
OE1001PB	1000
OE2001PB	2000
OE3001PB	3000
OE5001PB	5000



EVA INFUSION BAGS w. NEEDLE - FREE CONNECTOR

Ethylene vinyl acetate bags containing a needle-free connector at the port end. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.



TRANSPARENT EVA BAGS w. NEEDLE-FREE CONNECTOR	
REF NO	VOLUME (mL)
OE250V	250
OE500V	500
OE1000V	1000
OE2000V	2000
OE3000V	3000
OE5000V	5000


UV PROTECTED EVA BAGS w. NEEDLE-FREE CONNECTOR	
REF NO	VOLUME (mL)
OE250PV	250
OE500PV	500
OE1000PV	1000
OE2000PV	2000
OE3000PV	3000
OE5000PV	5000



EVA INFUSION BAGS w. EXTENSION LINE

Ethylene vinyl acetate bags containing an extension line. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.



TRANSPARENT EVA BAGS w. EXTENSION LINE & FLL CONNECTOR	
REF NO	VOLUME (mL)
OE250-11	250
OE500-11	500
OE1000-11	1000
OE2000-11	2000
OE3000-11	3000
OE5000-11	5000

TRANSPARENT EVA BAGS w. EXTENSION LINE & BIG BORE CONNECTOR	
REF NO	VOLUME (mL)
OE250B-11	250
OE500B-11	500
OE1000B-11	1000
OE2000B-11	2000
OE3000B-11	3000
OE5000B-11	5000

UV PROTECTED EVA BAGS w. EXTENSION LINE & FLL CONNECTOR	
REF NO	VOLUME (mL)
OE250P-12	250
OE500P-12	500
OE1000P-12	1000
OE2000P-12	2000
OE3000P-12	3000
OE5000P-12	5000


UV PROTECTED EVA BAGS w. EXTENSION LINE & BIG BORE CONNECTOR	
REF NO	VOLUME (mL)
OE250PB-12	250
OE500PB-12	500
OE1000PB-12	1000
OE2000PB-12	2000
OE3000PB-12	3000
OE5000PB-12	5000



EVA INFUSION BAGS w. 3 - WAY TRANSFER SET

Ethylene vinyl acetate bags designed for separate storage of lipid emulsion or other special solutions. UV protected options are available. Available in different volumes from 250mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 250mL - 500mL - 1000mL - 2000mL - 3000mL - 5000mL volume options.





EVA INFUSION BAGS w. 3 - WAY TRANSFER SET

TRANSPARENT EVA BAGS w. 3-WAY TRANSFER SET & FLL CONNECTOR

REF NO	VOLUME (mL)
OE250-30	250
OE500-30	500
OE1000-30	1000
OE2000-30	2000
OE3000-30	3000
OE5000-30	5000

UV PROTECTED EVA BAGS w. 3-WAY TRANSFER SET & FLL CONNECTOR

REF NO	VOLUME (mL)
OE250P-30	250
OE500P-30	500
OE1000P-30	1000
OE2000P-30	2000
OE3000P-30	3000
OE5000P-30	5000

TRANSPARENT EVA BAGS w. 3-WAY TRANSFER SET & MLL CONNECTOR

REF NO	VOLUME (mL)
OE251-30	250
OE501-30	500
OE1001-30	1000
OE2001-30	2000
OE3001-30	3000
OE5001-30	5000

UV PROTECTED EVA BAGS w. 3-WAY TRANSFER SET & MLL CONNECTOR

REF NO	VOLUME (mL)
OE251P-30	250
OE501P-30	500
OE1001P-30	1000
OE2001P-30	2000
OE3001P-30	3000
OE5001P-30	5000

TRANSPARENT EVA BAGS w. 3-WAY TRANSFER SET & BIG BORE CONNECTOR

REF NO	VOLUME (mL)
OE250B-30	250
OE500B-30	500
OE1000B-30	1000
OE2000B-30	2000
OE3000B-30	3000
OE5000B-30	5000

UV PROTECTED EVA BAGS w. 3-WAY TRANSFER SET & BIG BORE CONNECTOR

REF NO	VOLUME (mL)
OE250PB-30	250
OE500PB-30	500
OE1000PB-30	1000
OE2000PB-30	2000
OE3000PB-30	3000
OE5000PB-30	5000

TRANSPARENT EVA BAGS w. 3-WAY TRANSFER SET & NEEDLE-FREE CONNECTOR

REF NO	VOLUME (mL)
OE250V-30	250
OE500V-30	500
OE1000V-30	1000
OE2000V-30	2000
OE3000V-30	3000
OE5000V-30	5000

UV PROTECTED EVA BAGS w. 3-WAY TRANSFER SET & NEEDLE-FREE CONNECTOR

REF NO	VOLUME (mL)
OE250PV-30	250
OE500PV-30	500
OE1000PV-30	1000
OE2000PV-30	2000
OE3000PV-30	3000
OE5000PV-30	5000



EVA INFUSION BAGS w. 6 - WAY TRANSFER SET

Ethylene vinyl acetate bags designed for separated storage of lipid emulsion or other special solutions. UV protected options are available. Available in different volumes from 1000mL to 5000mL. Bag clamps are designed as openable or non-openable on demand.

- DEHP and LATEX free.
- Suitable to use with automatic compounders.
- UV protected options.
- 1.000 mL - 2.000 mL - 3.000 mL
5.000 mL volume options.



TRANSPARENT EVA BAGS w. 6-WAY TRANSFER SET & FLL CONNECTOR

REF NO	VOLUME (mL)
OE1000-40	1000
OE2000-40	2000
OE3000-40	3000
OE5000-40	5000

TRANSPARENT EVA BAGS w. 6-WAY TRANSFER SET & MLL CONNECTOR

REF NO	VOLUME (mL)
OE1001-40	1000
OE2001-40	2000
OE3001-40	3000
OE5001-40	5000

TRANSPARENT EVA BAGS w. 6-WAY TRANSFER SET & BIG BORE CONNECTOR

REF NO	VOLUME (mL)
OE1000B-40	1000
OE2000B-40	2000
OE3000B-40	3000
OE5000B-40	5000

IUV PROTECTED EVA BAGS w. 6-WAY TRANSFER SET & FLL CONNECTOR

REF NO	VOLUME (mL)
OE1000P-40	1000
OE2000P-40	2000
OE3000P-40	3000
OE5000P-40	5000

IUV PROTECTED EVA BAGS w. 6-WAY TRANSFER SET & MLL CONNECTOR

REF NO	VOLUME (mL)
OE1001P-40	1000
OE2001P-40	2000
OE3001P-40	3000
OE5001P-40	5000


IUV PROTECTED EVA BAGS w. 6-WAY TRANSFER SET & BIG BORE CONNECTOR

REF NO	VOLUME (mL)
OE1000PB-40	1000
OE2000PB-40	2000
OE3000PB-40	3000
OE5000PB-40	5000



NUTRITION

SETS



NUTRITION

SETS

ENTERAL NUTRITION SETS – 38

ENTERAL ADMINISTRATION SETS – 40

ENTERAL EXTENSION SETS – 41

NEONATAL PARENTERAL NUTRITION SETS – 42

PEDIATRIC SETS w. 0,2 MICRON FILTER – 43

PEDIATRIC SETS w. 1,2 MICRON FILTER – 44

PEDIATRIC SETS w. 1,2 MICRON FILTER - 160 CM – 45


TWO - WAY PEDIATRIC SETS w. 0,2 MICRON FILTER – 46

TWO - WAY PEDIATRIC SETS w. 1,2 MICRON FILTER – 47

ADULT PARENTERAL NUTRITION SETS – 48

ADULT SETS w. 0,2 MICRON FILTER – 49


ADULT SETS w. 1,2 MICRON FILTER – 49



ENTERAL NUTRITION SETS

NUTRIERA Enteral Nutrition Sets are designed for feeding patients who cannot take nutrients orally. It has a transparent bag structure where the nutritional liquid can be seen easily, an administration line for transferring the liquid, and a label area where patient information can be written.

- DEHP and LATEX free.
- Includes Enfit connector end with a cap.
- Includes luer lock connector compatible with PEG connectors/jejunal probes.
- 500mL - 800mL - 1000mL - 1200mL - 1500mL - 1800mL - 2000mL bag options.



ENTERAL NUTRITION SETS



ENTERAL NUTRITION SETS

SETS USED w. FEEDING PUMP	VOLUME (mL)
EN210P1 Enteral Nutrition Set	500
EN220P1 Enteral Nutrition Set	800
EN230P1 Enteral Nutrition Set	1000
EN240P1 Enteral Nutrition Set	1200
EN250P1 Enteral Nutrition Set	1500
EN260P1 Enteral Nutrition Set	1800
EN270P1 Enteral Nutrition Set	2000

ENTERAL NUTRITION SETS - DOUBLE BAG


SETS USED w. FEEDING PUMP	VOLUME (mL)
EN234P1 Enteral Nutrition Set w. Flush Bag	1000
EN244P1 Enteral Nutrition Set w. Flush Bag	1200
EN254P1 Enteral Nutrition Set w. Flush Bag	1500

ENTERAL NUTRITION SETS

SETS USED w. GRAVITY	VOLUME (mL)
EN210G1 Enteral Nutrition Set	500
EN220G1 Enteral Nutrition Set	800
EN230G1 Enteral Nutrition Set	1000
EN240G1 Enteral Nutrition Set	1200
EN250G1 Enteral Nutrition Set	1500
EN260G1 Enteral Nutrition Set	1800
EN270G1 Enteral Nutrition Set	2000

ENTERAL NUTRITION SETS - BOTTLE CAPPED

ENTERAL NUTRITION SETS - BOTTLE CAPPED
EN185P1 Enteral Nutrition Set w. Bottle Cap



ENTERAL ADMINISTRATION SETS

NUTRIERA Enteral Administration Sets are designed for feeding patients who cannot take nutrients orally. Sets are used by connecting to enteral bags. Contains an administration line for liquid transfer and a label where the patient information will be written.

- DEHP and LATEX free.
- Includes Enfit connector end with a cap.
- Includes luer lock connector compatible with PEG connectors/jejunal probes.
- Available to use with enteral feeding pump or gravity flow.

SETS USED w. FEEDING PUMP

EN180P1 Enteral Administration Set

SETS USED w. GRAVITY


EN180G1 Enteral Administration Set




ENTERAL EXTENSION SETS

NUTRIERA Enteral Extension Sets are available in varying lengths according to user needs. It can be connected to the Y connector of enteral nutrition sets and enteral administration sets.

- DEHP and LATEX free.
- Includes Enfit connector end with a cap.
- Includes luer lock connector compatible with PEG connectors/jejunal probes.




ENTERAL EXTENSION SETS	LENGTH (cm)
EN301U1 Enteral Extension Set	15
EN401U1 Enteral Extension Set	30
EN501U1 Enteral Extension Set	60
EN601U1 Enteral Extension Set	90
EN701U1 Enteral Extension Set	100
EN801U1 Enteral Extension Set	120




NEONATAL PARENTERAL NUTRITION SETS

NUTRIERA Neonatal Parenteral Nutrition Sets are configured in accordance with the principles of TPN administration in pediatrics. Contains 0,2 μm or 1,2 μm filter. Allows for proper filtering of nutritional content and lipid content separately to support two-way administration.

- Provides removal of air bubbles in liquid infusion.
- Removes large particles, bacteria and endotoxins from the system.
- Prevents harmful effects of daylight with its UV light-protected structure.
- Supports patient safety.



NEONATAL PARENTERAL NUTRITION SETS PEDIATRIC SETS w. 0,2 MICRON FILTER




PVC SETS (DEHP - FREE)


REF NO	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV100P1.102	✓	-	1,5X2,5	✓
IV100S1.102	✓	-	1,5X2,5	-
IV101P1.102	✓	-	1,5X2,7	✓
IV101S1.102	✓	-	1,5X2,7	-
IV102P1.102	✓	-	1,5X3,0	✓
IV102S1.102	✓	-	1,5X3,0	-
IV110P1.102	✓	✓	1,5X2,5	✓
IV110S1.102	✓	✓	1,5X2,5	-
IV111P1.102	✓	✓	1,5X2,7	✓
IV111S1.102	✓	✓	1,5X2,7	-
IV112P1.102	✓	✓	1,5X3,0	✓
IV112S1.102	✓	✓	1,5X3,0	-

POLYURETHANE SETS

REF NO	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV100P0.102	✓	-	1,5X2,5	✓
IV100S0.102	✓	-	1,5X2,5	-
IV101P0.102	✓	-	1,5X2,7	✓
IV101S0.102	✓	-	1,5X2,7	-
IV102P0.102	✓	-	1,5X3,0	✓
IV102S0.102	✓	-	1,5X3,0	-
IV110P0.102	✓	✓	1,5X2,5	✓
IV110S0.102	✓	✓	1,5X2,5	-
IV111P0.102	✓	✓	1,5X2,7	✓
IV111S0.102	✓	✓	1,5X2,7	-
IV112P0.102	✓	✓	1,5X3,0	✓
IV112S0.102	✓	✓	1,5X3,0	-



NEONATAL PARENTERAL NUTRITION SETS PEDIATRIC SETS w. 1,2 MICRON FILTER




PVC SETS (DEHP - FREE)

REF NO	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV100P1.112	✓	-	1,5X2,5	✓
IV100S1.112	✓	-	1,5X2,5	-
IV101P1.112	✓	-	1,5X2,7	✓
IV101S1.112	✓	-	1,5X2,7	-
IV102P1.112	✓	-	1,5X3,0	✓
IV102S1.112	✓	-	1,5X3,0	-
IV110P1.112	✓	✓	1,5X2,5	✓
IV110S1.112	✓	✓	1,5X2,5	-
IV111P1.112	✓	✓	1,5X2,7	✓
IV111S1.112	✓	✓	1,5X2,7	-
IV112P1.112	✓	✓	1,5X3,0	✓
IV112S1.112	✓	✓	1,5X3,0	-

POLYURETHANE SETS

REF NO	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV100P0.112	✓	-	1,5X2,5	✓
IV100S0.112	✓	-	1,5X2,5	-
IV101P0.112	✓	-	1,5X2,7	✓
IV101S0.112	✓	-	1,5X2,7	-
IV102P0.112	✓	-	1,5X3,0	✓
IV102S0.112	✓	-	1,5X3,0	-
IV110P0.112	✓	✓	0,65X2,5	✓
IV110S0.112	✓	✓	0,65X2,5	-
IV111P0.112	✓	✓	1,5X2,7	✓
IV111S0.112	✓	✓	1,5X2,7	-
IV112P0.112	✓	✓	1,5X3,0	✓
IV112S0.112	✓	✓	1,5X3,0	-

NEONATAL PARENTERAL NUTRITION SETS PEDIATRIC SETS w. 1,2 MICRON FILTER - 160 CM




PVC SETS (DEHP - FREE)

REF NO	1,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV100P1.112-160	✓	1,5X2,5	✓
IV100S1.112-160	✓	1,5X2,5	-
IV101P1.112-160	✓	1,5X2,7	✓
IV101S1.112-160	✓	1,5X2,7	-
IV102P1.112-160	✓	1,5X3,0	✓
IV102S1.112-160	✓	1,5X3,0	-
IV104P1.112-160	✓	0,65X2,5	✓
IV104S1.112-160	✓	0,65X2,5	-


POLYURETHANE SETS

REF NO	1,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV100P0.112-160	✓	1,5X2,5	✓
IV100S0.112-160	✓	1,5X2,5	-
IV101P0.112-160	✓	1,5X2,7	✓
IV101S0.112-160	✓	1,5X2,7	-
IV102P0.112-160	✓	1,5X3,0	✓
IV102S0.112-160	✓	1,5X3,0	-
IV104P0.112-160	✓	0,65X2,5	✓
IV104S0.112-160	✓	0,65X2,5	-



NEONATAL PARENTERAL NUTRITION SETS

TWO - WAY PEDIATRIC SETS w. 0,2 MICRON FILTER




PVC SETS (DEHP - FREE)

REF NO	0,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV400P1.102	✓	1,5X2,5	✓
IV400S1.102	✓	1,5X2,5	-
IV401P1.102	✓	1,5X2,7	✓
IV401S1.102	✓	1,5X2,7	-
IV402P1.102	✓	1,5X3,0	✓
IV402S1.102	✓	1,5X3,0	-

POLYURETHANE SETS

REF NO	0,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV400P0.102	✓	1,5X2,5	✓
IV400S0.102	✓	1,5X2,5	-
IV401P0.102	✓	1,5X2,7	✓
IV401S0.102	✓	1,5X2,7	-
IV402P0.102	✓	1,5X3,0	✓
IV402S0.102	✓	1,5X3,0	-

NEONATAL PARENTERAL NUTRITION SETS TWO - WAY PEDIATRIC SETS w. 1,2 MICRON FILTER




PVC SETS (DEHP - FREE)

REF NO	1,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV400P1.112	✓	1,5X2,5	✓
IV400S1.112	✓	1,5X2,5	-
IV401P1.112	✓	1,5X2,7	✓
IV401S1.112	✓	1,5X2,7	-
IV402P1.112	✓	1,5X3,0	✓
IV402S1.112	✓	1,5X3,0	-

POLYURETHANE SETS


REF NO	1,2 MICRON FILTER	DIAMETER	UV PROTECTION
IV400P0.112	✓	1,5X2,5	✓
IV400S0.112	✓	1,5X2,5	-
IV401P0.112	✓	1,5X2,7	✓
IV401S0.112	✓	1,5X2,7	-
IV402P0.112	✓	1,5X3,0	✓
IV402S0.112	✓	1,5X3,0	-



ADULT PARENTERAL NUTRITION SETS

NUTRIERA Adult Parenteral Nutrition Sets are configured in accordance with the principles of TPN administration in adults. Contains 0,2 μm or 1,2 μm filter. Allows for proper filtering of nutritional content and lipid content separately to support two-way administration.

- Provides removal of air bubbles in liquid infusion.
- Removes large particles, bacteria and endotoxins from the system.
- Prevents the harmful effects of daylight with its UV light-protected structure.
- Supports patient safety.



ADULT PARENTERAL NUTRITION SETS ADULT SETS w. 0,2 MICRON FILTER



PVC SETS (DEHP - FREE)

REF NO	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV103P1.002	✓	-	3,0X4,1	✓
IV103S1.002	✓	-	3,0X4,1	-
IV113P1.002	✓	-	3,0X4,1	✓
IV113S1.002	✓	-	3,0X4,1	-

POLYURETHANE SETS

REF NO	0,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV103P0.002	✓	-	3,0X4,1	✓
IV103S0.002	✓	-	3,0X4,1	-
IV113P0.002	✓	-	3,0X4,1	✓
IV113S0.002	✓	-	3,0X4,1	-

ADULT PARENTERAL NUTRITION SETS ADULT SETS w. 1,2 MICRON FILTER



PVC SETS (DEHP - FREE)

REF NO	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV103P1.012	✓	-	3,0X4,1	✓
IV103S1.012	✓	-	3,0X4,1	-
IV113P1.012	✓	-	3,0X4,1	✓
IV113S1.012	✓	-	3,0X4,1	-

POLYURETHANE SETS

REF NO	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV103P0.012	✓	-	3,0X4,1	✓
IV103S0.012	✓	-	3,0X4,1	-
IV113P0.012	✓	-	3,0X4,1	✓
IV113S0.012	✓	-	3,0X4,1	-



INFUSION

THERAPY



INFUSION

INFUSION SETS – 52

EXTENSION SETS – 54

STANDARD EXTENSION SETS – **55**

EXTENSION SETS w. FILTER & NEEDLE-FREE CONNECTORS – **63**

EXTENSION SETS - PRESSURE RESISTANT TO 250 PSI – **65**

EXTENSION SETS - PRESSURE RESISTANT TO 600 PSI – **67**

EXTENSION SETS - PRESSURE RESISTANT TO 800 PSI – **68**

EXTENSION SETS - PRESSURE RESISTANT TO 1200 PSI – **70**

MANIFOLDS – **71**

EXTENSION SETS w. MANIFOLD – **72**

EXTENSION SETS w. MANIFOLD & VALVE – **73**

EXTENSION SETS w. 3-WAY STOPCOCK – **74**

INFUSION

THERAPY

NEEDLE - FREE CONNECTORS (VALVES) – 75

ONE - WAY VALVE SETS - 10 CM – **76**

2 - WAY VALVE SETS - 10 CM – **76**

2 - WAY VALVE SETS - 15 CM – **76**

3 - WAY VALVE SETS - 10 CM – **77**

3 - WAY VALVE SETS - 15 CM – **77**

INFUSION ACCESSORIES – 78

HEPARIN CAP – **78**

TAMPER EVIDENT CAP – **78**

COMBI STOPPER – **78**

BAG SPIKE w. CHECK VALVE & NEEDLE-FREE CONNECTOR – **78**


BAG SPIKE w. NEEDLE-FREE CONNECTOR – **78**

LIQUID TRANSFER SET – **78**

TRANSFER SPIKE – **78**



INFUSION SETS



ONCOERA Infusion Sets are used for safe intravenous drug and fluid transfer to the patients in all wards, intensive care units and emergency departments that provide critical care where IV drug administration is carried out meticulously.

- DEHP and LATEX free.
- UV protected set options are available.
- Removes particles with a 15 µm particle filter in the drip chamber.
- Needle/needle-free Y port options are available.
- Special production can be made for appropriate infusion pumps.



LN001P1



LN002S1



LN007P1



LN007S1



LN077S1



INFUSION

INFUSION SETS

PVC SETS (DEHP - FREE)

REF NO	DRIP CHAMBER w. SPIKE	ROLLER CLAMP	Y INJECTION PORT	NEEDLE - FREE CONNECTOR	CLAMP	MALE LUER	UV PROTECTION	LENGTH
LN001S1	✓	✓	-	-	-	✓	-	230CM
LN002P1*	✓	✓	-	-	-	✓	✓	230CM
LN003P1	✓	✓	-	-	-	-	✓	195CM
LN006P1	✓	✓	✓	-	-	✓	-	200CM
LN007P1	✓	✓	✓	-	✓	✓	✓	225CM
LN007S1	✓	✓	✓	-	✓	✓	-	230CM
LN015P1	✓	✓	✓	-	-	-	✓	180CM
LN077P1	✓	✓	✓	-	✓	✓	✓	225CM
LN077S1	✓	✓	✓	-	✓	✓	-	225CM
LN107S1	-	✓	-	-	✓	✓	✓	180CM

*Contains silicone line. Silicone line optimizes the sensitivity of the administration set to deliver fluid through the infusion pump. Gravity infusion sets do not have a silicone line.

POLYURETHANE SETS

REF NO	DRIP CHAMBER w. SPIKE	ROLLER CLAMP	Y INJECTION PORT	NEEDLE - FREE CONNECTOR	CLAMP	MALE LUER	UV PROTECTION	LENGTH
LN001S2	✓	✓	-	-	-	✓	-	230CM
LN002P2*	✓	✓	-	-	-	✓	✓	230CM
LN003P2	✓	✓	-	-	-	-	✓	195CM
LN006P2	✓	✓	✓	-	-	✓	-	200CM
LN007P2	✓	✓	✓	-	✓	✓	✓	225CM
LN007S2	✓	✓	✓	-	✓	✓	-	230CM
LN015P2	✓	✓	✓	-	-	-	✓	180CM
LN077P2	✓	✓	✓	-	✓	✓	✓	225CM
LN077S2	✓	✓	✓	-	✓	✓	-	225CM
LN107S2	-	✓	-	-	✓	✓	✓	180CM


*Contains silicone line. Silicone line optimizes the sensitivity of the administration set to deliver fluid through the infusion pump. Gravity infusion sets do not have a silicone line.



EXTENSION SETS

ONCOERA Extension Sets are frequently used for the administration of infusion solutions and medications. The sets expand existing infusion systems for specific requirements of administration as needed and can also be used with other medical devices such as syringe pumps. ONCOERA offers a wide range of extension sets. The sets are optionally equipped with an integrated check valve and female/male luer lock adapter. Extension sets are available for a variety of lengths, diameters and materials.

- DEHP and LATEX free.
- It is pressure resistant.
- UV protected set options are available.
- Options with filter & needle-free connector are available.





EXTENSION SETS STANDARD EXTENSION SETS MLL/FLL Extension Sets

UV PROTECTED PVC SETS (DEHP - FREE)

LENGTH - CM	MLL / FLL 1x2mm	MLL / FLL 1x2,5mm	MLL / FLL 1,5x2,5mm	MLL / FLL 1,5x2,7mm	MLL / FLL 1,5x3mm	MLL / FLL 2,5x4,1mm	MLL / FLL 3x4,1mm
15	LN621PO-15A	LN621PO-15B	LN621PO-15C	LN621PO-15D	LN621PO-15E	LN621PO-15K	LN621PO-15L
30	LN621PO-30A	LN621PO-30B	LN621PO-30C	LN621PO-30D	LN621PO-30E	LN621PO-30K	LN621PO-30L
45	LN621PO-45A	LN621PO-45B	LN621PO-45C	LN621PO-45D	LN621PO-45E	LN621PO-45K	LN621PO-45L
60	LN621PO-60A	LN621PO-60B	LN621PO-60C	LN621PO-60D	LN621PO-60E	LN621PO-60K	LN621PO-60L
75	LN621PO-75A	LN621PO-75B	LN621PO-75C	LN621PO-75D	LN621PO-75E	LN621PO-75K	LN621PO-75L
90	LN621PO-90A	LN621PO-90B	LN621PO-90C	LN621PO-90D	LN621PO-90E	LN621PO-90K	LN621PO-90L
100	LN621PO-100A	LN621PO-100B	LN621PO-100C	LN621PO-100D	LN621PO-100E	LN621PO-100K	LN621PO-100L
120	LN621PO-120A	LN621PO-120B	LN621PO-120C	LN621PO-120D	LN621PO-120E	LN621PO-120K	LN621PO-120L
145	LN621PO-145A	LN621PO-145B	LN621PO-145C	LN621PO-145D	LN621PO-145E	LN621PO-145K	LN621PO-145L
150	LN621PO-150A	LN621PO-150B	LN621PO-150C	LN621PO-150D	LN621PO-150E	LN621PO-150K	LN621PO-150L
180	LN621PO-180A	LN621PO-180B	LN621PO-180C	LN621PO-180D	LN621PO-180E	LN621PO-180K	LN621PO-180L
200	LN621PO-200A	LN621PO-200B	LN621PO-200C	LN621PO-200D	LN621PO-200E	LN621PO-200K	LN621PO-200L
300	LN621PO-300A	LN621PO-300B	LN621PO-300C	LN621PO-300D	LN621PO-300E	LN621PO-300K	LN621PO-300L

TRANSPARENT PVC SETS (DEHP - FREE)

LENGTH - CM	MLL / FLL 1x2mm	MLL / FLL 1x2,5mm	MLL / FLL 1,5x2,5 mm	MLL / FLL 1,5x2,7mm	MLL / FLL 1,5x3mm	MLL / FLL 2,5x4,1mm	MLL / FLL 3x4,1mm
15	LN621SO-15A	LN621SO-15B	LN621SO-15C	LN621SO-15D	LN621SO-15E	LN621SO-15K	LN621SO-15L
30	LN621SO-30A	LN621SO-30B	LN621SO-30C	LN621SO-30D	LN621SO-30E	LN621SO-30K	LN621SO-30L
45	LN621SO-45A	LN621SO-45B	LN621SO-45C	LN621SO-45D	LN621SO-45E	LN621SO-45K	LN621SO-45L
60	LN621SO-60A	LN621SO-60B	LN621SO-60C	LN621SO-60D	LN621SO-60E	LN621SO-60K	LN621SO-60L
75	LN621SO-75A	LN621SO-75B	LN621SO-75C	LN621SO-75D	LN621SO-75E	LN621SO-75K	LN621SO-75L
90	LN621SO-90A	LN621SO-90B	LN621SO-90C	LN621SO-90D	LN621SO-90E	LN621SO-90K	LN621SO-90L
100	LN621SO-100A	LN621SO-100B	LN621SO-100C	LN621SO-100D	LN621SO-100E	LN621SO-100K	LN621SO-100L
120	LN621SO-120A	LN621SO-120B	LN621SO-120C	LN621SO-120D	LN621SO-120E	LN621SO-120K	LN621SO-120L
145	LN621SO-145A	LN621SO-145B	LN621SO-145C	LN621SO-145D	LN621SO-145E	LN621SO-145K	LN621SO-145L
150	LN621SO-150A	LN621SO-150B	LN621SO-150C	LN621SO-150D	LN621SO-150E	LN621SO-150K	LN621SO-150L
180	LN621SO-180A	LN621SO-180B	LN621SO-180C	LN621SO-180D	LN621SO-180E	LN621SO-180K	LN621SO-180L
200	LN621SO-200A	LN621SO-200B	LN621SO-200C	LN621SO-200D	LN621SO-200E	LN621SO-200K	LN621SO-200L
300	LN621SO-300A	LN621SO-300B	LN621SO-300C	LN621SO-300D	LN621SO-300E	LN621SO-300K	LN621SO-300L



EXTENSION SETS

STANDARD EXTENSION SETS

MLL/FLL Extension Sets

UV PROTECTED POLYURETHANE SETS

LENGTH - CM	MLL / FLL 1x2mm	MLL / FLL 1x2,5mm	MLL / FLL 1,5x2,5mm	MLL / FLL 1,5x2,7mm	MLL / FLL 1,5x3mm	MLL / FLL 2,5x4,1mm	MLL / FLL 3x4,1mm
15	LN721P0-15A	LN721P0-15B	LN721P0-15C	LN721P0-15D	LN721P0-15E	LN721P0-15K	LN721P0-15L
30	LN721P0-30A	LN721P0-30B	LN721P0-30C	LN721P0-30D	LN721P0-30E	LN721P0-30K	LN721P0-30L
45	LN721P0-45A	LN721P0-45B	LN721P0-45C	LN721P0-45D	LN721P0-45E	LN721P0-45K	LN721P0-45L
60	LN721P0-60A	LN721P0-60B	LN721P0-60C	LN721P0-60D	LN721P0-60E	LN721P0-60K	LN721P0-60L
75	LN721P0-75A	LN721P0-75B	LN721P0-75C	LN721P0-75D	LN721P0-75E	LN721P0-75K	LN721P0-75L
90	LN721P0-90A	LN721P0-90B	LN721P0-90C	LN721P0-90D	LN721P0-90E	LN721P0-90K	LN721P0-90L
100	LN721P0-100A	LN721P0-100B	LN721P0-100C	LN721P0-100D	LN721P0-100E	LN721P0-100K	LN721P0-100L
120	LN721P0-120A	LN721P0-120B	LN721P0-120C	LN721P0-120D	LN721P0-120E	LN721P0-120K	LN721P0-120L
145	LN721P0-145A	LN721P0-145B	LN721P0-145C	LN721P0-145D	LN721P0-145E	LN721P0-145K	LN721P0-145L
150	LN721P0-150A	LN721P0-150B	LN721P0-150C	LN721P0-150D	LN721P0-150E	LN721P0-150K	LN721P0-150L
180	LN721P0-180A	LN721P0-180B	LN721P0-180C	LN721P0-180D	LN721P0-180E	LN721P0-180K	LN721P0-180L
200	LN721P0-200A	LN721P0-200B	LN721P0-200C	LN721P0-200D	LN721P0-200E	LN721P0-200K	LN721P0-200L
300	LN721P0-300A	LN721P0-300B	LN721P0-300C	LN721P0-300D	LN721P0-300E	LN721P0-300K	LN721P0-300L

TRANSPARENT POLYURETHANE SETS

LENGTH - CM	MLL / FLL 1x2mm	MLL / FLL 1x2,5mm	MLL / FLL 1,5x2,5mm	MLL / FLL 1,5x2,7mm	MLL / FLL 1,5x3mm	MLL / FLL 2,5x4,1mm	MLL / FLL 3x4,1mm
15	LN721S0-15A	LN721S0-15B	LN721S0-15C	LN721S0-15D	LN721S0-15E	LN721S0-15K	LN721S0-15L
30	LN721S0-30A	LN721S0-30B	LN721S0-30C	LN721S0-30D	LN721S0-30E	LN721S0-30K	LN721S0-30L
45	LN721S0-45A	LN721S0-45B	LN721S0-45C	LN721S0-45D	LN721S0-45E	LN721S0-45K	LN721S0-45L
60	LN721S0-60A	LN721S0-60B	LN721S0-60C	LN721S0-60D	LN721S0-60E	LN721S0-60K	LN721S0-60L
75	LN721S0-75A	LN721S0-75B	LN721S0-75C	LN721S0-75D	LN721S0-75E	LN721S0-75K	LN721S0-75L
90	LN721S0-90A	LN721S0-90B	LN721S0-90C	LN721S0-90D	LN721S0-90E	LN721S0-90K	LN721S0-90L
100	LN721S0-100A	LN721S0-100B	LN721S0-100C	LN721S0-100D	LN721S0-100E	LN721S0-100K	LN721S0-100L
120	LN721S0-120A	LN721S0-120B	LN721S0-120C	LN721S0-120D	LN721S0-120E	LN721S0-120K	LN721S0-120L
145	LN721S0-145A	LN721S0-145B	LN721S0-145C	LN721S0-145D	LN721S0-145E	LN721S0-145K	LN721S0-145L
150	LN721S0-150A	LN721S0-150B	LN721S0-150C	LN721S0-150D	LN721S0-150E	LN721S0-150K	LN721S0-150L
180	LN721S0-180A	LN721S0-180B	LN721S0-180C	LN721S0-180D	LN721S0-180E	LN721S0-180K	LN721S0-180L
200	LN721S0-200A	LN721S0-200B	LN721S0-200C	LN721S0-200D	LN721S0-200E	LN721S0-200K	LN721S0-200L
300	LN721S0-300A	LN721S0-300B	LN721S0-300C	LN721S0-300D	LN721S0-300E	LN721S0-300K	LN721S0-300L



EXTENSION SETS

STANDARD EXTENSION SETS

MLL/MLL Extension Sets

UV PROTECTED PVC SETS (DEHP - FREE)

LENGTH - CM	MLL / MLL 1x2mm	MLL / MLL 1x2,5mm	MLL / MLL 1,5x2,5mm	MLL / MLL 1,5x2,7mm	MLL / MLL 1,5x3mm	MLL / MLL 2,5x4,1mm	MLL / MLL 3x4,1mm
15	LN622PO-15A	LN622PO-15B	LN622PO-15C	LN622PO-15D	LN622PO-15E	LN622PO-15K	LN622PO-15L
30	LN622PO-30A	LN622PO-30B	LN622PO-30C	LN622PO-30D	LN622PO-30E	LN622PO-30K	LN622PO-30L
45	LN622PO-45A	LN622PO-45B	LN622PO-45C	LN622PO-45D	LN622PO-45E	LN622PO-45K	LN622PO-45L
60	LN622PO-60A	LN622PO-60B	LN622PO-60C	LN622PO-60D	LN622PO-60E	LN622PO-60K	LN622PO-60L
75	LN622PO-75A	LN622PO-75B	LN622PO-75C	LN622PO-75D	LN622PO-75E	LN622PO-75K	LN622PO-75L
90	LN622PO-90A	LN622PO-90B	LN622PO-90C	LN622PO-90D	LN622PO-90E	LN622PO-90K	LN622PO-90L
100	LN622PO-100A	LN622PO-100B	LN622PO-100C	LN622PO-100D	LN622PO-100E	LN622PO-100K	LN622PO-100L
120	LN622PO-120A	LN622PO-120B	LN622PO-120C	LN622PO-120D	LN622PO-120E	LN622PO-120K	LN622PO-120L
145	LN622PO-145A	LN622PO-145B	LN622PO-145C	LN622PO-145D	LN622PO-145E	LN622PO-145K	LN622PO-145L
150	LN622PO-150A	LN622PO-150B	LN622PO-150C	LN622PO-150D	LN622PO-150E	LN622PO-150K	LN622PO-150L
180	LN622PO-180A	LN622PO-180B	LN622PO-180C	LN622PO-180D	LN622PO-180E	LN622PO-180K	LN622PO-180L
200	LN622PO-200A	LN622PO-200B	LN622PO-200C	LN622PO-200D	LN622PO-200E	LN622PO-200K	LN622PO-200L
300	LN622PO-300A	LN622PO-300B	LN622PO-300C	LN622PO-300D	LN622PO-300E	LN622PO-300K	LN622PO-300L

TRANSPARENT PVC SETS (DEHP - FREE)

LENGTH - CM	MLL / MLL 1x2mm	MLL / MLL 1x2,5mm	MLL / MLL 1,5x2,5mm	MLL / MLL 1,5x2,7mm	MLL / MLL 1,5x3mm	MLL / MLL 2,5x4,1mm	MLL / MLL 3x4,1mm
15	LN622SO-15A	LN622SO-15B	LN622SO-15C	LN622SO-15D	LN622SO-15E	LN622SO-15K	LN622SO-15L
30	LN622SO-30A	LN622SO-30B	LN622SO-30C	LN622SO-30D	LN622SO-30E	LN622SO-30K	LN622SO-30L
45	LN622SO-45A	LN622SO-45B	LN622SO-45C	LN622SO-45D	LN622SO-45E	LN622SO-45K	LN622SO-45L
60	LN622SO-60A	LN622SO-60B	LN622SO-60C	LN622SO-60D	LN622SO-60E	LN622SO-60K	LN622SO-60L
75	LN622SO-75A	LN622SO-75B	LN622SO-75C	LN622SO-75D	LN622SO-75E	LN622SO-75K	LN622SO-75L
90	LN622SO-90A	LN622SO-90B	LN622SO-90C	LN622SO-90D	LN622SO-90E	LN622SO-90K	LN622SO-90L
100	LN622SO-100A	LN622SO-100B	LN622SO-100C	LN622SO-100D	LN622SO-100E	LN622SO-100K	LN622SO-100L
120	LN622SO-120A	LN622SO-120B	LN622SO-120C	LN622SO-120D	LN622SO-120E	LN622SO-120K	LN622SO-120L
145	LN622SO-145A	LN622SO-145B	LN622SO-145C	LN622SO-145D	LN622SO-145E	LN622SO-145K	LN622SO-145L
150	LN622SO-150A	LN622SO-150B	LN622SO-150C	LN622SO-150D	LN622SO-150E	LN622SO-150K	LN622SO-150L
180	LN622SO-180A	LN622SO-180B	LN622SO-180C	LN622SO-180D	LN622SO-180E	LN622SO-180K	LN622SO-180L
200	LN622SO-200A	LN622SO-200B	LN622SO-200C	LN622SO-200D	LN622SO-200E	LN622SO-200K	LN622SO-200L
300	LN622SO-300A	LN622SO-300B	LN622SO-300C	LN622SO-300D	LN622SO-300E	LN622SO-300K	LN622SO-300L



EXTENSION SETS

STANDARD EXTENSION SETS

MLL/MLL Extension Sets

UV PROTECTED POLYURETHANE SETS

LENGTH - CM	MLL / MLL 1x2mm	MLL / MLL 1x2,5mm	MLL / MLL 1,5x2,5mm	MLL / MLL 1,5x2,7mm	MLL / MLL 1,5x3mm	MLL / MLL 2,5x4,1mm	MLL / MLL 3x4,1mm
15	LN722PO-15A	LN722PO-15B	LN722PO-15C	LN722PO-15D	LN722PO-15E	LN722PO-15K	LN722PO-15L
30	LN722PO-30A	LN722PO-30B	LN722PO-30C	LN722PO-30D	LN722PO-30E	LN722PO-30K	LN722PO-30L
45	LN722PO-45A	LN722PO-45B	LN722PO-45C	LN722PO-45D	LN722PO-45E	LN722PO-45K	LN722PO-45L
60	LN722PO-60A	LN722PO-60B	LN722PO-60C	LN722PO-60D	LN722PO-60E	LN722PO-60K	LN722PO-60L
75	LN722PO-75A	LN722PO-75B	LN722PO-75C	LN722PO-75D	LN722PO-75E	LN722PO-75K	LN722PO-75L
90	LN722PO-90A	LN722PO-90B	LN722PO-90C	LN722PO-90D	LN722PO-90E	LN722PO-90K	LN722PO-90L
100	LN722PO-100A	LN722PO-100B	LN722PO-100C	LN722PO-100D	LN722PO-100E	LN722PO-100K	LN722PO-100L
120	LN722PO-120A	LN722PO-120B	LN722PO-120C	LN722PO-120D	LN722PO-120E	LN722PO-120K	LN722PO-120L
145	LN722PO-145A	LN722PO-145B	LN722PO-145C	LN722PO-145D	LN722PO-145E	LN722PO-145K	LN722PO-145L
150	LN722PO-150A	LN722PO-150B	LN722PO-150C	LN722PO-150D	LN722PO-150E	LN722PO-150K	LN722PO-150L
180	LN722PO-180A	LN722PO-180B	LN722PO-180C	LN722PO-180D	LN722PO-180E	LN722PO-180K	LN722PO-180L
200	LN722PO-200A	LN722PO-200B	LN722PO-200C	LN722PO-200D	LN722PO-200E	LN722PO-200K	LN722PO-200L
300	LN722PS0-300A	LN722PO-300B	LN722PO-300C	LN722PO-300D	LN722PO-300E	LN722PO-300K	LN722PO-300L

TRANSPARENT POLYURETHANE SETS

LENGTH - CM	MLL / MLL 1x2mm	MLL / MLL 1x2,5mm	MLL / MLL 1,5x2,5mm	MLL / MLL 1,5x2,7mm	MLL / MLL 1,5x3mm	MLL / MLL 2,5x4,1mm	MLL / MLL 3x4,1mm
15	LN722S0-15A	LN722S0-15B	LN722S0-15C	LN722S0-15D	LN722S0-15E	LN722S0-15K	LN722S0-15L
30	LN722S0-30A	LN722S0-30B	LN722S0-30C	LN722S0-30D	LN722S0-30E	LN722S0-30K	LN722S0-30L
45	LN722S0-45A	LN722S0-45B	LN722S0-45C	LN722S0-45D	LN722S0-45E	LN722S0-45K	LN722S0-45L
60	LN722S0-60A	LN722S0-60B	LN722S0-60C	LN722S0-60D	LN722S0-60E	LN722S060K	LN722S0-60L
75	LN722S0-75A	LN722S0-75B	LN722S0-75C	LN722S0-75D	LN722S0-75E	LN722S0-75K	LN722S0-75L
90	LN722S0-90A	LN722S0-90B	LN722S0-90C	LN722S0-90D	LN722S0-90E	LN722S0-90K	LN722S0-90L
100	LN722S0-100A	LN722S0-100B	LN722S0-100C	LN722S0-100D	LN722S0-100E	LN722S0-100K	LN722S0-100L
120	LN722S0-120A	LN722S0-120B	LN722S0-120C	LN722S0-120D	LN722S0-120E	LN722S0-120K	LN722S0-120L
145	LN722S0-145A	LN722S0-145B	LN722S0-145C	LN722S0-145D	LN722S0-145E	LN722S0-145K	LN722S0-145L
150	LN722S0-150A	LN722S0-150B	LN722S0-150C	LN722S0-150D	LN722S0-150E	LN722S0-150K	LN722S0-150L
180	LN722S0-180A	LN722S0-180B	LN722S0-180C	LN722S0-180D	LN722S0-180E	LN722S0-180K	LN722S0-180L
200	LN722S0-200A	LN722S0-200B	LN722S0-200C	LN722S0-200D	LN722S0-200E	LN722S0-200K	LN722S0-200L
300	LN722S0-300A	LN722S0-300B	LN722S0-300C	LN722S0-300D	LN722S0-300E	LN722S0-300K	LN722S0-300L



EXTENSION SETS STANDARD EXTENSION SETS FLL/FLL Extension Sets

UV PROTECTED PVC SETS (DEHP - FREE)

LENGTH - CM	FLL / FLL 1x2mm	FLL / FLL 1x2,5mm	FLL / FLL 1,5x2,5mm	FLL / FLL 1,5x2,7mm	FLL / FLL 1,5x3mm	FLL / FLL 2,5x4,1mm	FLL / FLL 3x4,1mm
15	LN623PO-15A	LN623PO-15B	LN623PO-15C	LN623PO-15D	LN623PO-15E	LN623PO-15K	LN623PO-15L
30	LN623PO-30A	LN623PO-30B	LN623PO-30C	LN623PO-30D	LN623PO-30E	LN623PO-30K	LN623PO-30L
45	LN623PO-45A	LN623PO-45B	LN623PO-45C	LN623PO-45D	LN623PO-45E	LN623PO-45K	LN623PO-45L
60	LN623PO-60A	LN623PO-60B	LN623PO-60C	LN623PO-60D	LN623PO-60E	LN623PO-60K	LN623PO-60L
75	LN623PO-75A	LN623PO-75B	LN623PO-75C	LN623PO-75D	LN623PO-75E	LN623PO-75K	LN623PO-75L
90	LN623PO-90A	LN623PO-90B	LN623PO-90C	LN623PO-90D	LN623PO-90E	LN623PO-90K	LN623PO-90L
100	LN623PO-100A	LN623PO-100B	LN623PO-100C	LN623PO-100D	LN623PO-100E	LN623PO-100K	LN623PO-100L
120	LN623PO-120A	LN623PO-120B	LN623PO-120C	LN623PO-120D	LN623PO-120E	LN623PO-120K	LN623PO-120L
145	LN623PO-145A	LN623PO-145B	LN623PO-145C	LN623PO-145D	LN623PO-145E	LN623PO-145K	LN623PO-145L
150	LN623PO-150A	LN623PO-150B	LN623PO-150C	LN623PO-150D	LN623PO-150E	LN623PO-150K	LN623PO-150L
180	LN623PO-180A	LN623PO-180B	LN623PO-180C	LN623PO-180D	LN623PO-180E	LN623PO-180K	LN623PO-180L
200	LN623PO-200A	LN623PO-200B	LN623PO-200C	LN623PO-200D	LN623PO-200E	LN623PO-200K	LN623PO-200L
300	LN623PO-300A	LN623PO-300B	LN623PO-300C	LN623PO-300D	LN623PO-300E	LN623PO-300K	LN623PO-300L

TRANSPARENT PVC SETS (DEHP - FREE)

LENGTH - CM	FLL / FLL 1x2mm	FLL / FLL 1x2,5mm	FLL / FLL 1,5x2,5mm	FLL / FLL 1,5x2,7mm	FLL / FLL 1,5x3mm	FLL / FLL 2,5x4,1mm	FLL / FLL 3x4,1mm
15	LN623SO-15A	LN623SO-15B	LN623SO-15C	LN623SO-15D	LN623SO-15E	LN623SO-15K	LN623SO-15L
30	LN623SO-30A	LN623SO-30B	LN623SO-30C	LN623SO-30D	LN623SO-30E	LN623SO-30K	LN623SO-30L
45	LN623SO-45A	LN623SO-45B	LN623SO-45C	LN623SO-45D	LN623SO-45E	LN623SO-45K	LN623SO-45L
60	LN623SO-60A	LN623SO-60B	LN623SO-60C	LN623SO-60D	LN623SO-60E	LN623SO-60K	LN623SO-60L
75	LN623SO-75A	LN623SO-75B	LN623SO-75C	LN623SO-75D	LN623SO-75E	LN623SO-75K	LN623SO-75L
90	LN623SO-90A	LN623SO-90B	LN623SO-90C	LN623SO-90D	LN623SO-90E	LN623SO-90K	LN623SO-90L
100	LN623SO-100A	LN623SO-100B	LN623SO-100C	LN623SO-100D	LN623SO-100E	LN623SO-100K	LN623SO-100L
120	LN623SO-120A	LN623SO-120B	LN623SO-120C	LN623SO-120D	LN623SO-120E	LN623SO-120K	LN623SO-120L
145	LN623SO-145A	LN623SO-145B	LN623SO-145C	LN623SO-145D	LN623SO-145E	LN623SO-145K	LN623SO-145L
150	LN623SO-150A	LN623SO-150B	LN623SO-150C	LN623SO-150D	LN623SO-150E	LN623SO-150K	LN623SO-150L
180	LN623SO-180A	LN623SO-180B	LN623SO-180C	LN623SO-180D	LN623SO-180E	LN623SO-180K	LN623SO-180L
200	LN623SO-200A	LN623SO-200B	LN623SO-200C	LN623SO-200D	LN623SO-200E	LN623SO-200K	LN623SO-200L
300	LN623SO-300A	LN623SO-300B	LN623SO-300C	LN623SO-300D	LN623SO-300E	LN623SO-300K	LN623SO-300L



EXTENSION SETS

STANDARD EXTENSION SETS

FLL/FLL Extension Sets

UV PROTECTED POLYURETHANE SETS

LENGTH - CM	FLL / FLL 1x2mm	FLL / FLL 1x2,5mm	FLL / FLL 1,5x2,5mm	FLL / FLL 1,5x2,7mm	FLL / FLL 1,5x3mm	FLL / FLL 2,5x4,1mm	FLL / FLL 3x4,1mm
15	LN723PO-15A	LN723PO-15B	LN723PO-15C	LN723PO-15D	LN723PO-15E	LN723PO-15K	LN723PO-15L
30	LN723PO-30A	LN723PO-30B	LN723PO-30C	LN723PO-30D	LN723PO-30E	LN723PO-30K	LN723PO-30L
45	LN723PO-45A	LN723PO-45B	LN723PO-45C	LN723PO-45D	LN723PO-45E	LN723PO-45K	LN723PO-45L
60	LN723PO-60A	LN723PO-60B	LN723PO-60C	LN723PO-60D	LN723PO-60E	LN723PO60K	LN723PO-60L
75	LN723PO-75A	LN723PO-75B	LN723PO-75C	LN723PO-75D	LN723PO-75E	LN723PO-75K	LN723PO-75L
90	LN723PO-90A	LN723PO-90B	LN723PO-90C	LN723PO-90D	LN723PO-90E	LN723PO-90K	LN723PO-90L
100	LN723PO-100A	LN723PO-100B	LN723PO-100C	LN723PO-100D	LN723PO-100E	LN723PO-100K	LN723PO-100L
120	LN723PO-120A	LN723PO-120B	LN723PO-120C	LN723PO-120D	LN723PO-120E	LN723PO-120K	LN723PO-120L
145	LN723PO-145A	LN723PO-145B	LN723PO-145C	LN723PO-145D	LN723PO-145E	LN723PO-145K	LN723PO-145L
150	LN723PO-150A	LN723PO-150B	LN723PO-150C	LN723PO-150D	LN723PO-150E	LN722S3PO-150K	LN723PO-150L
180	LN723PO-180A	LN723PO-180B	LN723PO-180C	LN723PO-180D	LN723PO-180E	LN723PO-180K	LN723PO-180L
200	LN723PO-200A	LN723PO-200B	LN723PO-200C	LN723PO-200D	LN723PO-200E	LN723PO-200K	LN723PO-200L
300	LN723PO-300A	LN723PO-300B	LN723PO-300C	LN723PO-300D	LN723PO-300E	LN723PO-300K	LN723PO-300L

TRANSPARENT POLYURETHANE SETS

LENGTH - CM	FLL / FLL 1x2mm	FLL / FLL 1x2,5mm	FLL / FLL 1,5x2,5mm	FLL / FLL 1,5x2,7mm	FLL / FLL 1,5x3mm	FLL / FLL 2,5x4,1mm	FLL / FLL 3x4,1mm
15	LN723SO-15A	LN723SO-15B	LN723SO-15C	LN723SO-15D	LN723SO-15E	LN723SO-15K	LN723SO-15L
30	LN723SO-30A	LN723SO-30B	LN723SO-30C	LN723SO-30D	LN723SO-30E	LN723SO-30K	LN723SO-30L
45	LN723SO-45A	LN723SO-45B	LN723SO-45C	LN723SO-45D	LN723SO-45E	LN723SO-45K	LN723SO-45L
60	LN723SO-60A	LN723SO-60B	LN723SO-60C	LN723SO-60D	LN723SO-60E	LN723SO60K	LN723SO-60L
75	LN723SO-75A	LN723SO-75B	LN723SO-75C	LN723SO-75D	LN723SO-75E	LN723SO-75K	LN723SO-75L
90	LN723SO-90A	LN723SO-90B	LN723SO-90C	LN723SO-90D	LN723SO-90E	LN723SO-90K	LN723SO-90L
100	LN723SO-100A	LN723SO-100B	LN723SO-100C	LN723SO-100D	LN723SO-100E	LN723SO-100K	LN723SO-100L
120	LN723SO-120A	LN723SO-120B	LN723SO-120C	LN723SO-120D	LN723SO-120E	LN723SO-120K	LN723SO-120L
145	LN723SO-145A	LN723SO-145B	LN723SO-145C	LN723SO-145D	LN723SO-145E	LN723SO-145K	LN723SO-145L
150	LN723SO-150A	LN723SO-150B	LN723SO-150C	LN723SO-150D	LN723SO-150E	LN722S3SO-150K	LN723SO-150L
180	LN723SO-180A	LN723SO-180B	LN723SO-180C	LN723SO-180D	LN723SO-180E	LN723SO-180K	LN723SO-180L
200	LN723SO-200A	LN723SO-200B	LN723SO-200C	LN723SO-200D	LN723SO-200E	LN723SO-200K	LN723SO-200L
300	LN723SO-300A	LN723SO-300B	LN723SO-300C	LN723SO-300D	LN723SO-300E	LN723SO-300K	LN723SO-300L



EXTENSION SETS STANDARD EXTENSION SETS ROTATING MLL/FLL Extension Sets

UV PROTECTED PVC SETS (DEHP - FREE)

LENGTH - CM	ROT MLL/FLL 1x2mm	ROT MLL/FLL 1x2,5mm	ROT MLL/FLL 1,5x2,5mm	ROT MLL/FLL 1,5x2,7mm	ROT MLL/FLL 1,5x3mm	ROT MLL/FLL 2,5X4,1mm	ROT MLL/FLL 3x4,1mm
15	LN624PO-15A	LN624PO-15B	LN624PO-15C	LN624PO-15D	LN624PO-15E	LN624PO-15K	LN624PO-15L
30	LN624PO-30A	LN624PO-30B	LN624PO-30C	LN624PO-30D	LN624PO-30E	LN624PO-30K	LN624PO-30L
45	LN624PO-45A	LN624PO-45B	LN624PO-45C	LN624PO-45D	LN624PO-45E	LN624PO-45K	LN624PO-45L
60	LN624PO-60A	LN624PO-60B	LN624PO-60C	LN624PO-60D	LN624PO-60E	LN624PO-60K	LN624PO-60L
75	LN624PO-75A	LN624PO-75B	LN624PO-75C	LN624PO-75D	LN624PO-75E	LN624PO-75K	LN624PO-75L
90	LN624PO-90A	LN624PO-90B	LN624PO-90C	LN624PO-90D	LN624PO-90E	LN624PO-90K	LN624PO-90L
100	LN624PO-100A	LN624PO-100B	LN624PO-100C	LN624PO-100D	LN624PO-100E	LN624PO-100K	LN624PO-100L
120	LN624PO-120A	LN624PO-120B	LN624PO-120C	LN624PO-120D	LN624PO-120E	LN624PO-120K	LN624PO-120L
145	LN624PO-145A	LN624PO-145B	LN624PO-145C	LN624PO-145D	LN624PO-145E	LN624PO-145K	LN624PO-145L
150	LN624PO-150A	LN624PO-150B	LN624PO-150C	LN624PO-150D	LN624PO-150E	LN624PO-150K	LN624PO-150L
180	LN624PO-180A	LN624PO-180B	LN624PO-180C	LN624PO-180D	LN624PO-180E	LN624PO-180K	LN624PO-180L
200	LN624PO-200A	LN624PO-200B	LN624PO-200C	LN624PO-200D	LN624PO-200E	LN624PO-200K	LN624PO-200L
300	LN624PO-300A	LN624PO-300B	LN624PO-300C	LN624PO-300D	LN624PO-300E	LN624PO-300K	LN624PO-300L

TRANSPARENT PVC SETS (DEHP - FREE)

LENGTH - CM	ROT MLL/FLL 1x2mm	ROT MLL/FLL 1x2,5mm	ROT MLL/FLL 1,5x2,5mm	ROT MLL/FLL 1,5x2,7mm	ROT MLL/FLL 1,5x3mm	ROT MLL/FLL 2,5X4,1mm	ROT MLL/FLL 3x4,1mm
15	LN624SO-15A	LN624SO-15B	LN624SO-15C	LN624SO-15D	LN624SO-15E	LN624SO-15K	LN624SO-15L
30	LN624SO-30A	LN624SO-30B	LN624SO-30C	LN624SO-30D	LN624SO-30E	LN624SO-30K	LN624SO-30L
45	LN624SO-45A	LN624SO-45B	LN624SO-45C	LN624SO-45D	LN624SO-45E	LN624SO-45K	LN624SO-45L
60	LN624SO-60A	LN624SO-60B	LN624SO-60C	LN624SO-60D	LN624SO-60E	LN624SO-60K	LN624SO-60L
75	LN624SO-75A	LN624SO-75B	LN624SO-75C	LN624SO-75D	LN624SO-75E	LN624SO-75K	LN624SO-75L
90	LN624SO-90A	LN624SO-90B	LN624SO-90C	LN624SO-90D	LN624SO-90E	LN624SO-90K	LN624SO-90L
100	LN624SO-100A	LN624SO-100B	LN624SO-100C	LN624SO-100D	LN624SO-100E	LN624SO-100K	LN624SO-100L
120	LN624SO-120A	LN624SO-120B	LN624SO-120C	LN624SO-120D	LN624SO-120E	LN624SO-120K	LN624SO-120L
145	LN624SO-145A	LN624SO-145B	LN624SO-145C	LN624SO-145D	LN624SO-145E	LN624SO-145K	LN624SO-145L
150	LN624SO-150A	LN624SO-150B	LN624SO-150C	LN624SO-150D	LN624SO-150E	LN624SO-150K	LN624SO-150L
180	LN624SO-180A	LN624SO-180B	LN624SO-180C	LN624SO-180D	LN624SO-180E	LN624SO-180K	LN624SO-180L
200	LN624SO-200A	LN624SO-200B	LN624SO-200C	LN624SO-200D	LN624SO-200E	LN624SO-200K	LN624SO-200L
300	LN624SO-300A	LN624SO-300B	LN624SO-300C	LN624SO-300D	LN624SO-300E	LN624SO-300K	LN624SO-300L



EXTENSION SETS

STANDARD EXTENSION SETS

ROTATING MLL/FLL Extension Sets

UV PROTECTED POLYURETHANE SETS

LENGTH - CM	ROT MLL/FLL 1x2mm	ROT MLL/FLL 1x2,5mm	ROT MLL/FLL 1,5x2,5mm	ROT MLL/FLL 1,5x2,7mm	ROT MLL/FLL 1,5x3mm	ROT MLL/FLL 2,5X4,1mm	ROT MLL/FLL 3x4,1mm
15	LN724PO-15A	LN724PO-15B	LN724PO-15C	LN724PO-15D	LN724PO-15E	LN724PO-15K	LN724PO-15L
30	LN724PO-30A	LN724PO-30B	LN724PO-30C	LN724PO-30D	LN724PO-30E	LN724PO-30K	LN724PO-30L
45	LN724PO-45A	LN724PO-45B	LN724PO-45C	LN724PO-45D	LN724PO-45E	LN724PO-45K	LN724PO-45L
60	LN724PO-60A	LN724PO-60B	LN724PO-60C	LN724PO-60D	LN724PO-60E	LN724PO-60K	LN724PO-60L
75	LN724PO-75A	LN724PO-75B	LN724PO-75C	LN724PO-75D	LN724PO-75E	LN724PO-75K	LN724PO-75L
90	LN724PO-90A	LN724PO-90B	LN724PO-90C	LN724PO-90D	LN724PO-90E	LN724PO-90K	LN724PO-90L
100	LN724PO-100A	LN724PO-100B	LN724PO-100C	LN724PO-100D	LN724PO-100E	LN724PO-100K	LN724PO-100L
120	LN724PO-120A	LN724PO-120B	LN724PO-120C	LN724PO-120D	LN724PO-120E	LN724PO-120K	LN724PO-120L
145	LN724PO-145A	LN724PO-145B	LN724PO-145C	LN724PO-145D	LN724PO-145E	LN724PO-145K	LN724PO-145L
150	LN724PO-150A	LN724PO-150B	LN724PO-150C	LN724PO-150D	LN724PO-150E	LN724PO-150K	LN724PO-150L
180	LN724PO-180A	LN724PO-180B	LN724PO-180C	LN724PO-180D	LN724PO-180E	LN724PO-180K	LN724PO-180L
200	LN724PO-200A	LN724PO-200B	LN724PO-200C	LN724PO-200D	LN724PO-200E	LN724PO-200K	LN724PO-200L
300	LN724PO-300A	LN724PO-300B	LN724PO-300C	LN724PO-300D	LN724PO-300E	LN724PO-300K	LN724PO-300L

TRANSPARENT POLYURETHANE SETS

LENGTH - CM	ROT MLL/FLL 1x2mm	ROT MLL/FLL 1x2,5mm	ROT MLL/FLL 1,5x2,5mm	ROT MLL/FLL 1,5x2,7mm	ROT MLL/FLL 1,5x3mm	ROT MLL/FLL 2,5X4,1mm	ROT MLL/FLL 3x4,1mm
15	LN724SO-15A	LN724SO-15B	LN724SO-15C	LN724SO-15D	LN724SO-15E	LN724SO-15K	LN724SO-15L
30	LN724SO-30A	LN724SO-30B	LN724SO-30C	LN724SO-30D	LN724SO-30E	LN724SO-30K	LN724SO-30L
45	LN724SO-45A	LN724SO-45B	LN724SO-45C	LN724SO-45D	LN724SO-45E	LN724SO-45K	LN724SO-45L
60	LN724SO-60A	LN724SO-60B	LN724SO-60C	LN724SO-60D	LN724SO-60E	LN724SO-60K	LN724SO-60L
75	LN724SO-75A	LN724SO-75B	LN724SO-75C	LN724SO-75D	LN724SO-75E	LN724SO-75K	LN724SO-75L
90	LN724SO-90A	LN724SO-90B	LN724SO-90C	LN724SO-90D	LN724SO-90E	LN724SO-90K	LN724SO-90L
100	LN724SO-100A	LN724SO-100B	LN724SO-100C	LN724SO-100D	LN724SO-100E	LN724SO-100K	LN724SO-100L
120	LN724SO-120A	LN724SO-120B	LN724SO-120C	LN724SO-120D	LN724SO-120E	LN724SO-120K	LN724SO-120L
145	LN724SO-145A	LN724SO-145B	LN724SO-145C	LN724SO-145D	LN724SO-145E	LN724SO-145K	LN724SO-145L
150	LN724SO-150A	LN724SO-150B	LN724SO-150C	LN724SO-150D	LN724SO-150E	LN724SO-150K	LN724SO-150L
180	LN724SO-180A	LN724SO-180B	LN724SO-180C	LN724SO-180D	LN724SO-180E	LN724SO-180K	LN724SO-180L
200	LN724SO-200A	LN724SO-200B	LN724SO-200C	LN724SO-200D	LN724SO-200E	LN724SO-200K	LN724SO-200L
300	LN724SO-300A	LN724SO-300B	LN724SO-300C	LN724SO-300D	LN724SO-300E	LN724SO-300K	LN724SO-300L



EXTENSION SETS

EXTENSION SETS w. FILTER & NEEDLE - FREE CONNECTORS
Extension Sets w. Two Filters & Two Needle - Free Connectors



PVC SETS (DEHP - FREE)

REF NO	0,2 MICRON FILTER	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV210P1.102	2 PCS	-	2 PCS	1,5X2,5	✓
IV210S1.102	2 PCS	-	2 PCS	1,5X2,5	-
IV211P1.102	2 PCS	-	2 PCS	1,5X2,7	✓
IV211S1.102	2 PCS	-	2 PCS	1,5X2,7	-
IV212P1.102	2 PCS	-	2 PCS	1,5X3,0	✓
IV212S1.102	2 PCS	-	2 PCS	1,5X3,0	-
IV210P1.112	-	2 PCS	2 PCS	1,5X2,5	✓
IV210S1.112	-	2 PCS	2 PCS	1,5X2,5	-
IV211P1.112	-	2 PCS	2 PCS	1,5X2,7	✓
IV211S1.112	-	2 PCS	2 PCS	1,5X2,7	-
IV212P1.112	-	2 PCS	2 PCS	1,5X3,0	✓
IV212S1.112	-	2 PCS	2 PCS	1,5X3,0	-

POLYURETHANE SETS

REF NO	0,2 MICRON FILTER	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV210P0.102	2 PCS	-	2 PCS	1,5X2,5	✓
IV210S0.102	2 PCS	-	2 PCS	1,5X2,5	-
IV211P0.102	2 PCS	-	2 PCS	1,5X2,7	✓
IV211S0.102	2 PCS	-	2 PCS	1,5X2,7	-
IV212P0.102	2 PCS	-	2 PCS	1,5X3,0	✓
IV212S0.102	2 PCS	-	2 PCS	1,5X3,0	-
IV210P0.112	-	2 PCS	2 PCS	1,5X2,5	✓
IV210S0.112	-	2 PCS	2 PCS	1,5X2,5	-
IV211P0.112	-	2 PCS	2 PCS	1,5X2,7	✓
IV211S0.112	-	2 PCS	2 PCS	1,5X2,7	-
IV212P0.112	-	2 PCS	2 PCS	1,5X3,0	✓
IV212S0.112	-	2 PCS	2 PCS	1,5X3,0	-



EXTENSION SETS

EXTENSION SETS w. FILTER & NEEDLE - FREE CONNECTORS

Extension Sets w. a Filter & Two Needle - Free Connectors



PVC SETS (DEHP - FREE)

REF NO	0,2 MICRON FILTER	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV410P1.102	1 PC	-	2 PCS	1,5X2,5	✓
IV410S1.102	1 PC	-	2 PCS	1,5X2,5	-
IV411P1.102	1 PC	-	2 PCS	1,5X2,7	✓
IV411S1.102	1 PC	-	2 PCS	1,5X2,7	-
IV412P1.102	1 PC	-	2 PCS	1,5X3,0	✓
IV412S1.102	1 PC	-	2 PCS	1,5X3,0	-
IV410P1.112	-	1 PC	2 PCS	1,5X2,5	✓
IV410S1.112	-	1 PC	2 PCS	1,5X2,5	-
IV411P1.112	-	1 PC	2 PCS	1,5X2,7	✓
IV411S1.112	-	1 PC	2 PCS	1,5X2,7	-
IV412P1.112	-	1 PC	2 PCS	1,5X3,0	✓
IV412S1.112	-	1 PC	2 PCS	1,5X3,0	-

POLYURETHANE SETS

REF NO	0,2 MICRON FILTER	1,2 MICRON FILTER	NEEDLE - FREE CONNECTOR	DIAMETER	UV PROTECTION
IV410PO.102	1 PC	-	2 PCS	1,5X2,5	✓
IV410SO.102	1 PC	-	2 PCS	1,5X2,5	-
IV411PO.102	1 PC	-	2 PCS	1,5X2,7	✓
IV411SO.102	1 PC	-	2 PCS	1,5X2,7	-
IV412PO.102	1 PC	-	2 PCS	1,5X3,0	✓
IV412SO.102	1 PC	-	2 PCS	1,5X3,0	-
IV410PO.112	-	1 PC	2 PCS	1,5X2,5	✓
IV410SO.112	-	1 PC	2 PCS	1,5X2,5	-
IV411PO.112	-	1 PC	2 PCS	1,5X2,7	✓
IV411SO.112	-	1 PC	2 PCS	1,5X2,7	-
IV412PO.112	-	1 PC	2 PCS	1,5X3,0	✓
IV412SO.112	-	1 PC	2 PCS	1,5X3,0	-



EXTENSION SETS

EXTENSION SETS - PRESSURE RESISTANT TO 250 PSI

1,5x2,5mm UV Protected & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x2,5mm	MLL / MLL 1,5x2,5mm	FLL / FLL 1,5x2,5mm	ROT MLL / FLL 1,5x2,5mm
15	LN621P1-15C	LN622P1-15C	LN623P1-15C	LN624P1-15C
30	LN621P1-30C	LN622P1-30C	LN623P1-30C	LN624P1-30C
45	LN621P1-45C	LN622P1-45C	LN623P1-45C	LN624P1-45C
60	LN621P1-60C	LN622P1-60C	LN623P1-60C	LN624P1-60C
75	LN621P1-75C	LN622P1-75C	LN623P1-75C	LN624P1-75C
90	LN621P1-90C	LN622P1-90C	LN623P1-90C	LN624P1-90C
100	LN621P1-100C	LN622P1-100C	LN623P1-100C	LN624P1-100C
120	LN621P1-120C	LN622P1-120C	LN623P1-120C	LN624P1-120C
145	LN621P1-145C	LN622P1-145C	LN623P1-145C	LN624P1-145C
150	LN621P1-150C	LN622P1-150C	LN623P1-150C	LN624P1-150C
180	LN621P1-180C	LN622P1-180C	LN623P1-180C	LN624P1-180C
200	LN621P1-200C	LN622P1-200C	LN623P1-200C	LN624P1-200C
300	LN621P1-300C	LN622P1-300C	LN623P1-300C	LN624P1-300C

1,5x2,5mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x2,5mm	MLL / MLL 1,5x2,5mm	FLL / FLL 1,5x2,5mm	ROT MLL / FLL 1,5x2,5mm
15	LN621S1-15C	LN622S1-15C	LN623S1-15C	LN624S1-15C
30	LN621S1-30C	LN622S1-30C	LN623S1-30C	LN624S1-30C
45	LN621S1-45C	LN622S1-45C	LN623S1-45C	LN624S1-45C
60	LN621S1-60C	LN622S1-60C	LN623S1-60C	LN624S1-60C
75	LN621S1-75C	LN622S1-75C	LN623S1-75C	LN624S1-75C
90	LN621S1-90C	LN622S1-90C	LN623S1-90C	LN624S1-90C
100	LN621S1-100C	LN622S1-100C	LN623S1-100C	LN624S1-100C
120	LN621S1-120C	LN622S1-120C	LN623S1-120C	LN624S1-120C
145	LN621S1-145C	LN622S1-145C	LN623S1-145C	LN624S1-145C
150	LN621S1-150C	LN622S1-150C	LN623S1-150C	LN624S1-150C
180	LN621S1-180C	LN622S1-180C	LN623S1-180C	LN624S1-180C
200	LN621S1-200C	LN622S1-200C	LN623S1-200C	LN624S1-200C
300	LN621S1-300C	LN622S1-300C	LN623S1-300C	LN624S1-300C

1,5x2,7mm UV Protected & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x2,7mm	MLL / MLL 1,5x2,7mm	FLL / FLL 1,5x2,7mm	ROT MLL / FLL 1,5x2,7mm
15	LN621P1-15D	LN622P1-15D	LN623P1-15D	LN624P1-15D
30	LN621P1-30D	LN622P1-30D	LN623P1-30D	LN624P1-30D
45	LN621P1-45D	LN622P1-45D	LN623P1-45D	LN624P1-45D
60	LN621P1-60D	LN622P1-60D	LN623P1-60D	LN624P1-60D
75	LN621P1-75D	LN622P1-75D	LN623P1-75D	LN624P1-75D
90	LN621P1-90D	LN622P1-90D	LN623P1-90D	LN624P1-90D
100	LN621P1-100D	LN622P1-100D	LN623P1-100D	LN624P1-100D
120	LN621P1-120D	LN622P1-120D	LN623P1-120D	LN624P1-120D
145	LN621P1-145D	LN622P1-145D	LN623P1-145D	LN624P1-145D
150	LN621P1-150D	LN622P1-150D	LN623P1-150D	LN624P1-150D
180	LN621P1-180D	LN622P1-180D	LN623P1-180D	LN624P1-180D
200	LN621P1-200D	LN622P1-200D	LN623P1-200D	LN624P1-200D
300	LN621P1-300D	LN622P1-300D	LN623P1-300D	LN624P1-300D

1,5x2,7mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x2,7mm	MLL / MLL 1,5x2,7mm	FLL / FLL 1,5x2,7mm	ROT MLL / FLL 1,5x2,7mm
15	LN621S1-15D	LN622S1-15D	LN623S1-15D	LN624S1-15D
30	LN621S1-30D	LN622S1-30D	LN623S1-30D	LN624S1-30D
45	LN621S1-45D	LN622S1-45D	LN623S1-45D	LN624S1-45D
60	LN621S1-60D	LN622S1-60D	LN623S1-60D	LN624S1-60D
75	LN621S1-75D	LN622S1-75D	LN623S1-75D	LN624S1-75D
90	LN621S1-90D	LN622S1-90D	LN623S1-90D	LN624S1-90D
100	LN621S1-100D	LN622S1-100D	LN623S1-100D	LN624S1-100D
120	LN621S1-120D	LN622S1-120D	LN623S1-120D	LN624S1-120D
145	LN621S1-145D	LN622S1-145D	LN623S1-145D	LN624S1-145D
150	LN621S1-150D	LN622S1-150D	LN623S1-150D	LN624S1-150D
180	LN621S1-180D	LN622S1-180D	LN623S1-180D	LN624S1-180D
200	LN621S1-200D	LN622S1-200D	LN623S1-200D	LN624S1-200D
300	LN621S1-300D	LN622S1-300D	LN623S1-300D	LN624S1-300D



EXTENSION SETS

EXTENSION SETS - PRESSURE RESISTANT TO 250 PSI

1,5x3mm UV Protected & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x3mm	MLL / MLL 1,5x3mm	FLL / FLL 1,5x3mm	ROT MLL / FLL 1,5x3mm
15	LN621P1-15E	LN622P1-15E	LN623P1-15E	LN624P1-15E
30	LN621P1-30E	LN622P1-30E	LN623P1-30E	LN624P1-30E
45	LN621P1-45E	LN622P1-45E	LN623P1-45E	LN624P1-45E
60	LN621P1-60E	LN622P1-60E	LN623P1-60E	LN624P1-60E
75	LN621P1-75E	LN622P1-75E	LN623P1-75E	LN624P1-75E
90	LN621P1-90E	LN622P1-90E	LN623P1-90E	LN624P1-90E
100	LN621P1-100E	LN622P1-100E	LN623P1-100E	LN624P1-100E
120	LN621P1-120E	LN622P1-120E	LN623P1-120E	LN624P1-120E
145	LN621P1-145E	LN622P1-145E	LN623P1-145E	LN624P1-145E
150	LN621P1-150E	LN622P1-150E	LN623P1-150E	LN624P1-150E
180	LN621P1-180E	LN622P1-180E	LN623P1-180E	LN624P1-180E
200	LN621P1-200E	LN622P1-200E	LN623P1-200E	LN624P1-200E
300	LN621P1-300E	LN622P1-300E	LN623P1-300E	LN624P1-300E

1,5x3mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x3mm	MLL / MLL 1,5x3mm	FLL / FLL 1,5x3mm	ROT MLL / FLL 1,5x3mm
15	LN621S1-15E	LN622S1-15E	LN623S1-15E	LN624S1-15E
30	LN621S1-30E	LN622S1-30E	LN623S1-30E	LN624S1-30E
45	LN621S1-45E	LN622S1-45E	LN623S1-45E	LN624S1-45E
60	LN621S1-60E	LN622S1-60E	LN623S1-60E	LN624S1-60E
75	LN621S1-75E	LN622S1-75E	LN623S1-75E	LN624S1-75E
90	LN621S1-90E	LN622S1-90E	LN623S1-90E	LN624S1-90E
100	LN621S1-100E	LN622S1-100E	LN623S1-100E	LN624S1-100E
120	LN621S1-120E	LN622S1-120E	LN623S1-120E	LN624S1-120E
145	LN621S1-145E	LN622S1-145E	LN623S1-145E	LN624S1-145E
150	LN621S1-150E	LN622S1-150E	LN623S1-150E	LN624S1-150E
180	LN621S1-180E	LN622S1-180E	LN623S1-180E	LN624S1-180E
200	LN621S1-200E	LN622S1-200E	LN623S1-200E	LN624S1-200E
300	LN621S1-300E	LN622S1-300E	LN623S1-300E	LN624S1-300E

3x4,1mm UV Protected & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 3x4,1mm	MLL / MLL 3x4,1mm	FLL / FLL 3x4,1mm	ROT MLL / FLL 3x4,1mm
15	LN621P1-15L	LN622P1-15L	LN623P1-15L	LN624P1-15L
30	LN621P1-30L	LN622P1-30L	LN623P1-30L	LN624P1-30L
45	LN621P1-45L	LN622P1-45L	LN623P1-45L	LN624P1-45L
60	LN621P1-60L	LN622P1-60L	LN623P1-60L	LN624P1-60L
75	LN621P1-75L	LN622P1-75L	LN623P1-75L	LN624P1-75L
90	LN621P1-90L	LN622P1-90L	LN623P1-90L	LN624P1-90L
100	LN621P1-100L	LN622P1-100L	LN623P1-100L	LN624P1-100L
120	LN621P1-120L	LN622P1-120L	LN623P1-120L	LN624P1-120L
145	LN621P1-145L	LN622P1-145L	LN623P1-145L	LN624P1-145L
150	LN621P1-150L	LN622P1-150L	LN623P1-150L	LN624P1-150L
180	LN621P1-180L	LN622P1-180L	LN623P1-180L	LN624P1-180L
200	LN621P1-200L	LN622P1-200L	LN623P1-200L	LN624P1-200L
300	LN621P1-300L	LN622P1-300L	LN623P1-300L	LN624P1-300L

3x4,1mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 3x4,1mm	MLL / MLL 3x4,1mm	FLL / FLL 3x4,1mm	ROT MLL / FLL 3x4,1mm
15	LN621S1-15L	LN622S1-15L	LN623S1-15L	LN624S1-15L
30	LN621S1-30L	LN622S1-30L	LN623S1-30L	LN624S1-30L
45	LN621S1-45L	LN622S1-45L	LN623S1-45L	LN624S1-45L
60	LN621S1-60L	LN622S1-60L	LN623S1-60L	LN624S1-60L
75	LN621S1-75L	LN622S1-75L	LN623S1-75L	LN624S1-75L
90	LN621S1-90L	LN622S1-90L	LN623S1-90L	LN624S1-90L
100	LN621S1-100L	LN622S1-100L	LN623S1-100L	LN624S1-100L
120	LN621S1-120L	LN622S1-120L	LN623S1-120L	LN624S1-120L
145	LN621S1-145L	LN622S1-145L	LN623S1-145L	LN624S1-145L
150	LN621S1-150L	LN622S1-150L	LN623S1-150L	LN624S1-150L
180	LN621S1-180L	LN622S1-180L	LN623S1-180L	LN624S1-180L
200	LN621S1-200L	LN622S1-200L	LN623S1-200L	LN624S1-200L
300	LN621S1-300L	LN622S1-300L	LN623S1-300L	LN624S1-300L



EXTENSION SETS

EXTENSION SETS - PRESSURE RESISTANT TO 600 PSI

1,5x2,7mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x2,7mm	MLL / MLL 1,5x2,7mm	FLL / FLL 1,5x2,7mm	ROT MLL / FLL 1,5x2,7mm
15	LN631S2-15D	LN632S2-15D	LN633S2-15D	LN634S2-15D
30	LN631S2-30D	LN632S2-30D	LN633S2-30D	LN634S2-30D
45	LN631S2-45D	LN632S2-45D	LN633S2-45D	LN634S2-45D
60	LN631S2-60D	LN632S2-60D	LN633S2-60D	LN634S2-60D
75	LN631S2-75D	LN632S2-75D	LN633S2-75D	LN634S2-75D
90	LN631S2-90D	LN632S2-90D	LN633S2-90D	LN634S2-90D
100	LN631S2-100D	LN632S2-100D	LN633S2-100D	LN634S2-100D
120	LN631S2-120D	LN632S2-120D	LN633S2-120D	LN634S2-120D
145	LN631S2-145D	LN632S2-145D	LN633S2-145D	LN634S2-145D
150	LN631S2-150D	LN632S2-150D	LN633S2-150D	LN634S2-150D
180	LN631S2-180D	LN632S2-180D	LN633S2-180D	LN634S2-180D
200	LN631S2-200D	LN632S2-200D	LN633S2-200D	LN634S2-200D
300	LN631S2-300D	LN632S2-300D	LN633S2-300D	LN634S2-300D

1,5x3mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x3mm	MLL / MLL 1,5x3mm	FLL / FLL 1,5x3mm	ROT MLL / FLL 1,5x3mm
15	LN631S2-15E	LN632S2-15E	LN633S2-15E	LN634S2-15E
30	LN631S2-30E	LN632S2-30E	LN633S2-30E	LN634S2-30E
45	LN631S2-45E	LN632S2-45E	LN633S2-45E	LN634S2-45E
60	LN631S2-60E	LN632S2-60E	LN633S2-60E	LN634S2-60E
75	LN631S2-75E	LN632S2-75E	LN633S2-75E	LN634S2-75E
90	LN631S2-90E	LN632S2-90E	LN633S2-90E	LN634S2-90E
100	LN631S2-100E	LN632S2-100E	LN633S2-100E	LN634S2-100E
120	LN631S2-120E	LN632S2-120E	LN633S2-120E	LN634S2-120E
145	LN631S2-145E	LN632S2-145E	LN633S2-145E	LN634S2-145E
150	LN631S2-150E	LN632S2-150E	LN633S2-150E	LN634S2-150E
180	LN631S2-180E	LN632S2-180E	LN633S2-180E	LN634S2-180E
200	LN631S2-200E	LN632S2-200E	LN633S2-200E	LN634S2-200E
300	LN631S2-300E	LN632S2-300E	LN633S2-300E	LN634S2-300E

3x4,1mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 3X4,1mm	MLL / MLL 3X4,1mm	FLL / FLL 3X4,1mm	ROT MLL / FLL 3X4,1mm
15	LN631S2-15L	LN632S2-15L	LN633S2-15L	LN634S2-15L
30	LN631S2-30L	LN632S2-30L	LN633S2-30L	LN634S2-30L
45	LN631S2-45L	LN632S2-45L	LN633S2-45L	LN634S2-45L
60	LN631S2-60L	LN632S2-60L	LN633S2-60L	LN634S2-60L
75	LN631S2-75L	LN632S2-75L	LN633S2-75L	LN634S2-75L
90	LN631S2-90L	LN632S2-90L	LN633S2-90L	LN634S2-90L
100	LN631S2-100L	LN632S2-100L	LN633S2-100L	LN634S2-100L
120	LN631S2-120L	LN632S2-120L	LN633S2-120L	LN634S2-120L
145	LN631S2-145L	LN632S2-145L	LN633S2-145L	LN634S2-145L
150	LN631S2-150L	LN632S2-150L	LN633S2-150L	LN634S2-150L
180	LN631S2-180L	LN632S2-180L	LN633S2-180L	LN634S2-180L
200	LN631S2-200L	LN632S2-200L	LN633S2-200L	LN634S2-200L
300	LN631S2-300L	LN632S2-300L	LN633S2-300L	LN634S2-300L





EXTENSION SETS

EXTENSION SETS - PRESSURE RESISTANT TO 800 PSI

1,5x2,7mm UV Protected & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x2,7mm	MLL / MLL 1,5x2,7mm	FLL / FLL 1,5x2,7mm	ROT MLL / FLL 1,5x2,7mm
15	LN631P3-15D	LN632P3-15D	LN633P3-15D	LN634P3-15D
30	LN631P3-30D	LN632P3-30D	LN633P3-30D	LN634P3-30D
45	LN631P3-45D	LN632P3-45D	LN633P3-45D	LN634P3-45D
60	LN631P3-60D	LN632P3-60D	LN633P3-60D	LN634P3-60D
75	LN631P3-75D	LN632P3-75D	LN633P3-75D	LN634P3-75D
90	LN631P3-90D	LN632P3-90D	LN633P3-90D	LN634P3-90D
100	LN631P3-100D	LN632P3-100D	LN633P3-100D	LN634P3-100D
120	LN631P3-120D	LN632P3-120D	LN633P3-120D	LN634P3-120D
145	LN631P3-145D	LN632P3-145D	LN633P3-145D	LN634P3-145D
150	LN631P3-150D	LN632P3-150D	LN633P3-150D	LN634P3-150D
180	LN631P3-180D	LN632P3-180D	LN633P3-180D	LN634P3-180D
200	LN631P3-200D	LN632P3-200D	LN633P3-200D	LN634P3-200D
300	LN631P3-300D	LN632P3-300D	LN633P3-300D	LN634P3-300D

1,5x2,7mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x2,7mm	MLL / MLL 1,5x2,7mm	FLL / FLL 1,5x2,7mm	ROT MLL / FLL 1,5x2,7mm
15	LN631S3-15D	LN632S3-15D	LN633S3-15D	LN634S3-15D
30	LN631S3-30D	LN632S3-30D	LN633S3-30D	LN634S3-30D
45	LN631S3-45D	LN632S3-45D	LN633S3-45D	LN634S3-45D
60	LN631S3-60D	LN632S3-60D	LN633S3-60D	LN634S3-60D
75	LN631S3-75D	LN632S3-75D	LN633S3-75D	LN634S3-75D
90	LN631S3-90D	LN632S3-90D	LN633S3-90D	LN634S3-90D
100	LN631S3-100D	LN632S3-100D	LN633S3-100D	LN634S3-100D
120	LN631S3-120D	LN632S3-120D	LN633S3-120D	LN634S3-120D
145	LN631S3-145D	LN632S3-145D	LN633S3-145D	LN634S3-145D
150	LN631S3-150D	LN632S3-150D	LN633S3-150D	LN634S3-150D
180	LN631S3-180D	LN632S3-180D	LN633S3-180D	LN634S3-180D
200	LN631S3-200D	LN632S3-200D	LN633S3-200D	LN634S3-200D
300	LN631S3-300D	LN632S3-300D	LN633S3-300D	LN634S3-300D

1,5x3mm UV Protected & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x3mm	MLL / MLL 1,5x3mm	FLL / FLL 1,5x3mm	ROT MLL / FLL 1,5x3mm
15	LN631P3-15E	LN632P3-15E	LN633P3-15E	LN634P3-15E
30	LN631P3-30E	LN632P3-30E	LN633P3-30E	LN634P3-30E
45	LN631P3-45E	LN632P3-45E	LN633P3-45E	LN634P3-45E
60	LN631P3-60E	LN632P3-60E	LN633P3-60E	LN634P3-60E
75	LN631P3-75E	LN632P3-75E	LN633P3-75E	LN634P3-75E
90	LN631P3-90E	LN632P3-90E	LN633P3-90E	LN634P3-90E
100	LN631P3-100E	LN632P3-100E	LN633P3-100E	LN634P3-100E
120	LN631P3-120E	LN632P3-120E	LN633P3-120E	LN634P3-120E
145	LN631P3-145E	LN632P3-145E	LN633P3-145E	LN634P3-145E
150	LN631P3-150E	LN632P3-150E	LN633P3-150E	LN634P3-150E
180	LN631P3-180E	LN632P3-180E	LN633P3-180E	LN634P3-180E
200	LN631P3-200E	LN632P3-200E	LN633P3-200E	LN634P3-200E
300	LN631P3-300E	LN632P3-300E	LN633P3-300E	LN634P3-300E

1,5x3mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,5x3mm	MLL / MLL 1,5x3mm	FLL / FLL 1,5x3mm	ROT MLL / FLL 1,5x3mm
15	LN631S3-15E	LN632S3-15E	LN633S3-15E	LN634S3-15E
30	LN631S3-30E	LN632S3-30E	LN633S3-30E	LN634S3-30E
45	LN631S3-45E	LN632S3-45E	LN633S3-45E	LN634S3-45E
60	LN631S3-60E	LN632S3-60E	LN633S3-60E	LN634S3-60E
75	LN631S3-75E	LN632S3-75E	LN633S3-75E	LN634S3-75E
90	LN631S3-90E	LN632S3-90E	LN633S3-90E	LN634S3-90E
100	LN631S3-100E	LN632S3-100E	LN633S3-100E	LN634S3-100E
120	LN631S3-120E	LN632S3-120E	LN633S3-120E	LN634S3-120E
145	LN631S3-145E	LN632S3-145E	LN633S3-145E	LN634S3-145E
150	LN631S3-150E	LN632S3-150E	LN633S3-150E	LN634S3-150E
180	LN631S3-180E	LN632S3-180E	LN633S3-180E	LN634S3-180E
200	LN631S3-200E	LN632S3-200E	LN633S3-200E	LN634S3-200E
300	LN631S3-300E	LN632S3-300E	LN633S3-300E	LN634S3-300E



EXTENSION SETS

EXTENSION SETS - PRESSURE RESISTANT TO 800 PSI

2,3x4,1mm UV Protected & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 2,3X4,1mm	MLL / MLL 2,3X4,1mm	FLL / FLL 2,3X4,1mm	ROT MLL / FLL 2,3X4,1mm
15	LN631P3-15H	LN632P3-15H	LN633P3-15H	LN634P3-15H
30	LN631P3-30H	LN632P3-30H	LN633P3-30H	LN634P3-30H
45	LN631P3-45H	LN632P3-45H	LN633P3-45H	LN634P3-45H
60	LN631P3-60H	LN632P3-60H	LN633P3-60H	LN634P3-60H
75	LN631P3-75H	LN632P3-75H	LN633P3-75H	LN634P3-75H
90	LN631P3-90H	LN632P3-90H	LN633P3-90H	LN634P3-90H
100	LN631P3-100H	LN632P3-100H	LN633P3-100H	LN634P3-100H
120	LN631P3-120H	LN632P3-120H	LN633P3-120H	LN634P3-120H
145	LN631P3-145H	LN632P32-145H	LN633P3-145H	LN634P3-145H
150	LN631P3-150H	LN632P3-150H	LN633P3-150H	LN634P3-150H
180	LN631P3-180H	LN632P3-180H	LN633P3-180H	LN634P3-180H
200	LN631P3-200H	LN632P3-200H	LN633P3-200H	LN634P3-200H
300	LN631P3-300H	LN632P3-300H	LN633P3-300H	LN634P3-300H

2,3x4,1mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 2,3X4,1mm	MLL / MLL 2,3X4,1mm	FLL / FLL 2,3X4,1mm	ROT MLL / FLL 2,3X4,1mm
15	LN631S3-15H	LN632S3-15H	LN633S3-15H	LN634S3-15H
30	LN631S3-30H	LN632S3-30H	LN633S3-30H	LN634S3-30H
45	LN631S3-45H	LN632S3-45H	LN633S3-45H	LN634S3-45H
60	LN631S3-60H	LN632S3-60H	LN633S3-60H	LN634S3-60H
75	LN631S3-75H	LN632S3-75H	LN633S3-75H	LN634S3-75H
90	LN631S3-90H	LN632S3-90H	LN633S3-90H	LN634S3-90H
100	LN631S3-100H	LN632S3-100H	LN633S3-100H	LN634S3-100H
120	LN631S3-120H	LN632S3-120H	LN633S3-120H	LN634S3-120H
145	LN631S3-145H	LN632S32-145H	LN633S3-145H	LN634S3-145H
150	LN631S3-150H	LN632S3-150H	LN633S3-150H	LN634S3-150H
180	LN631S3-180H	LN632S3-180H	LN633S3-180H	LN634S3-180H
200	LN631S3-200H	LN632S3-200H	LN633S3-200H	LN634S3-200H
300	LN631S3-300H	LN632S3-300H	LN633S3-300H	LN634S3-300H

3x4,1mm UV Protected & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 3X4,1mm	MLL / MLL 3X4,1mm	FLL / FLL 3X4,1mm	ROT MLL / FLL 3X4,1mm
15	LN631P3-15L	LN632P3-15L	LN633P3-15L	LN634P3-15L
30	LN631P3-30L	LN632P3-30L	LN633P3-30L	LN634P3-30L
45	LN631P3-45L	LN632P3-45L	LN633P3-45L	LN634P3-45L
60	LN631P3-60L	LN632P3-60L	LN633P3-60L	LN634P3-60L
75	LN631P3-75L	LN632P3-75L	LN633P3-75L	LN634P3-75L
90	LN631P3-90L	LN632P3-90L	LN633P3-90L	LN634P3-90L
100	LN631P3-100L	LN632P3-100L	LN633P3-100L	LN634P3-100L
120	LN631P3-120L	LN632P3-120L	LN633P3-120L	LN634P3-120L
145	LN631P3-145L	LN632P32-145L	LN633P3-145L	LN634P3-145L
150	LN631P3-150L	LN632P3-150L	LN633P3-150L	LN634P3-150L
180	LN631P3-180L	LN632P3-180L	LN633P3-180L	LN634P3-180L
200	LN631P3-200L	LN632P3-200L	LN633P3-200L	LN634P3-200L
300	LN631P3-300L	LN632P3-300L	LN633P3-300L	LN634P3-300L


3x4,1mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 3X4,1mm	MLL / MLL 3X4,1mm	FLL / FLL 3X4,1mm	ROT MLL / FLL 3X4,1mm
15	LN631S3-15L	LN632S3-15L	LN633S3-15L	LN634S3-15L
30	LN631S3-30L	LN632S3-30L	LN633S3-30L	LN634S3-30L
45	LN631S3-45L	LN632S3-45L	LN633S3-45L	LN634S3-45L
60	LN631S3-60L	LN632S3-60L	LN633S3-60L	LN634S3-60L
75	LN631S3-75L	LN632S3-75L	LN633S3-75L	LN634S3-75L
90	LN631S3-90L	LN632S3-90L	LN633S3-90L	LN634S3-90L
100	LN631S3-100L	LN632S3-100L	LN633S3-100L	LN634S3-100L
120	LN631S3-120L	LN632S3-120L	LN633S3-120L	LN634S3-120L
145	LN631S3-145L	LN632S3-145L	LN633S3-145L	LN634S3-145L
150	LN631S3-150L	LN632S3-150L	LN633S3-150L	LN634S3-150L
180	LN631S3-180L	LN632S3-180L	LN633S3-180L	LN634S3-180L
200	LN631S3-200L	LN632S3-200L	LN633S3-200L	LN634S3-200L
300	LN631S3-300L	LN632S3-300L	LN633S3-300L	LN634S3-300L



EXTENSION SETS

EXTENSION SETS - PRESSURE RESISTANT TO 1200 PSI



1,6x4,1mm Transparent & PVC Extension Sets (DEHP - Free)


LENGTH (CM)	MLL / FLL 1,6X4,1mm	MLL / MLL 1,6X4,1mm	FLL / FLL 1,6X4,1mm	ROT MLL / FLL 1,6X4,1mm
15	LN631S4-15F	LN632S4-15F	LN633S4-15F	LN634S4-15F
30	LN631S4-30F	LN632S4-30F	LN633S4-30F	LN634S4-30F
45	LN631S4-45F	LN632S4-45F	LN633S4-45F	LN634S4-45F
60	LN631S4-60F	LN632S4-60F	LN633S4-60F	LN634S4-60F
75	LN631S4-75F	LN632S4-75F	LN633S4-75F	LN634S4-75F
90	LN631S4-90F	LN632S4-90F	LN633S4-90F	LN634S4-90F
100	LN631S4-100F	LN632S4-100F	LN633S4-100F	LN634S4-100F
120	LN631S4-120F	LN632S4-120F	LN633S4-120F	LN634S4-120F
145	LN631S4-145F	LN632S4-145F	LN633S4-145F	LN634S4-145F
150	LN631S4-150F	LN632S4-150F	LN633S4-150F	LN634S4-150F
180	LN631S4-180F	LN632S4-180F	LN633S4-180F	LN634S4-180F
200	LN631S4-200F	LN632S4-200F	LN633S4-200F	LN634S4-200F
300	LN631S4-300F	LN632S4-300F	LN633S4-300F	LN634S4-300F

1,7x4,9 Mm Transparent & PVC Extension Sets (DEHP - Free)

LENGTH (CM)	MLL / FLL 1,7X4,9mm	MLL / MLL 1,7X4,9mm	FLL / FLL 1,7X4,9mm	ROT MLL / FLL 1,7X4,9mm
15	LN631S4-15R	LN632S4-15R	LN633S4-15R	LN634S4-15R
30	LN631S4-30R	LN632S4-30R	LN633S4-30R	LN634S4-30R
45	LN631S4-45R	LN632S4-45R	LN633S4-45R	LN634S4-45R
60	LN631S4-60R	LN632S4-60R	LN633S4-60R	LN634S4-60R
75	LN631S4-75R	LN632S4-75R	LN633S4-75R	LN634S4-75R
90	LN631S4-90R	LN632S4-90R	LN633S4-90R	LN634S4-90R
100	LN631S4-100R	LN632S4-100R	LN633S4-100R	LN634S4-100R
120	LN631S4-120R	LN632S4-120R	LN633S4-120R	LN634S4-120R
145	LN631S4-145R	LN632S4-145R	LN633S4-145R	LN634S4-145R
150	LN631S4-150R	LN632S4-150R	LN633S4-150R	LN634S4-150R
180	LN631S4-180R	LN632S4-180R	LN633S4-180R	LN634S4-180R
200	LN631S4-200R	LN632S4-200R	LN633S4-200R	LN634S4-200R
300	LN631S4-300R	LN632S4-300R	LN633S4-300R	LN634S4-300R



EXTENSION SETS MANIFOLDS Manifold



MANIFOLD	
REF NO	DESCRIPTION
LN00100-002	2 - GANG MANIFOLD
LN00100-003	3 - GANG MANIFOLD
LN00100-004	4 - GANG MANIFOLD
LN00100-005	5 - GANG MANIFOLD

Manifold w. Valve




MANIFOLD w. VALVE	
REF NO	DESCRIPTION
LN00100-002V	2 - GANG MANIFOLD w. VALVE
LN00100-003V	3 - GANG MANIFOLD w. VALVE
LN00100-004V	4 - GANG MANIFOLD w. VALVE
LN00100-005V	5 - GANG MANIFOLD w. VALVE



EXTENSION SETS

EXTENSION SETS w. MANIFOLD



PVC SETS (DEHP - FREE)

LENGTH - CM	REF NO	DESCRIPTION
15	LN621S0-15K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
30	LN621S0-30K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
45	LN621S0-45K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
75	LN621S0-75K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
90	LN621S0-90K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
100	LN621S0-100K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
150	LN621S0-150K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
180	LN621S0-180K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
200	LN621S0-200K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
15	LN621S0-15K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
30	LN621S0-30K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
45	LN621S0-45K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
75	LN621S0-75K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
90	LN621S0-90K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
100	LN621S0-100K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
150	LN621S0-150K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
180	LN621S0-180K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
200	LN621S0-200K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5*4,1 mm


POLYURETHANE SETS

LENGTH - CM	REF NO	DESCRIPTION
15	LN721S0-15K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
30	LN721S0-30K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
45	LN721S0-45K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
75	LN721S0-75K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
90	LN721S0-90K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
100	LN721S0-100K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
150	LN721S0-150K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
180	LN721S0-180K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
200	LN721S0-200K3	EXTENSION SET w. 3 - GANG MANIFOLD 2,5x4,1 mm
15	LN721S0-15K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
30	LN721S0-30K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
45	LN721S0-45K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
75	LN721S0-75K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
90	LN721S0-90K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
100	LN721S0-100K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
150	LN721S0-150K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
180	LN721S0-180K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5x4,1 mm
200	LN721S0-200K5	EXTENSION SET w. 5 - GANG MANIFOLD 2,5*4,1 mm



EXTENSION SETS

EXTENSION SETS w. MANIFOLD & VALVE



PVC SETS (DEHP - FREE)

LENGTH - CM	REF NO	DESCRIPTION
15	LN621S0-15K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
30	LN621S0-30K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
45	LN621S0-45K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVED 2,5x4,1 mm
75	LN621S0-75K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
90	LN621S0-90K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
100	LN621S0-100K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
150	LN621S0-150K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
180	LN621S0-180K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
200	LN621S0-200K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
15	LN621S0-15K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
30	LN621S0-30K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
45	LN621S0-45K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
75	LN621S0-75K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVED 2,5x4,1 mm
90	LN621S0-90K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
100	LN621S0-100K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
150	LN621S0-150K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
180	LN621S0-180K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
200	LN621S0-200K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5*4,1 mm
300	LN621S0-300K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm

POLYURETHANE SETS

LENGTH - CM	REF NO	DESCRIPTION
15	LN721S0-15K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
30	LN721S0-30K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
45	LN721S0-45K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
75	LN721S0-75K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
90	LN721S0-90K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
100	LN721S0-100K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
150	LN721S0-150K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
180	LN721S0-180K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
200	LN721S0-200K3V	EXTENSION SET w. 3 - GANG MANIFOLD & VALVE 2,5x4,1 mm
15	LN721S0-15K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
30	LN721S0-30K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
45	LN721S0-45K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
75	LN721S0-75K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
90	LN721S0-90K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
100	LN721S0-100K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
150	LN721S0-150K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
180	LN721S0-180K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm
200	LN721S0-200K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5*4,1 mm
300	LN721S0-300K5V	EXTENSION SET w. 5 - GANG MANIFOLD & VALVE 2,5x4,1 mm



EXTENSION SETS

EXTENSION SETS w. 3 - WAY STOPCOCK



PVC SETS (DEHP - FREE)

LENGTH - CM	REF NO	DESCRIPTION
15	LN621S0-15K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
	LN621S0-30K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
45	LN621S0-45K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
	LN621S0-60K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
75	LN621S0-75K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
	LN621S0-90K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
100	LN621S0-100K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
	LN621S0-120K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
150	LN621S0-150K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
180	LN621S0-180K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
200	LN621S0-200K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
300	LN621S0-300K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm

POLYURETHANE SETS

LENGTH - CM	REF NO	DESCRIPTION
15	LN721S0-15K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
30	LN721S0-30K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
45	LN721S0-45K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
60	LN721S0-60K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
75	LN721S0-75K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
90	LN721S0-90K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
100	LN721S0-100K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
120	LN721S0-120K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
150	LN721S0-150K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
180	LN721S0-180K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
200	LN721S0-200K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm
300	LN721S0-300K1	EXTENSION SET w. 3 - WAY STOPCOCK 2,5x4,1mm



NEEDLE - FREE CONNECTORS (VALVES)



ONCOERA IV Needle-Free Connectors are the part of the split septum valve group. Needle-free connectors have become an integral part of intravenous therapy in recent years. The most important reason for this widespread acceptance is reduces the risk of infection of the patient and healthcare personnel. The CDC also recommends the use of split septum vein valves in intravenous therapy to reduce the risk of infection. Different extension set options are available, including single-valve, two-valve, three-valve and four-valve.

- Supports closed system.
- Protects compounded drug, personnel, and compounding environment.
- Suitable to use with all luer lock connectors.
- Forms a microbiological barrier for 7 days.
- Easy to disinfect.
- Prevents needle stick accidents.
- Provides low dead volume and high flow rate.
- Compatible with blood products and oncology drugs.
- DEHP - Free



IV201
VALVE w. TRANSPARENT SILICONE



IV301
VALVE w. BLUE SILICONE



NEEDLE - FREE CONNECTORS (VALVES) ONE - WAY VALVE SETS - 10 CM



1,5x2,5mm		1,5x2,7mm		1,5x3mm		3x4,1mm		1x2mm	
UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT
IV110P1	IV110S1	IV111P1	IV111S1	IV112P1	IV112S1	IV113P1	IV113S1	IV104P1.102	IV104S1.102
IV110P1-15	IV110S1-16	IV111P1-16	IV111S1-16	IV112P1-16	IV112S1-16	IV113P1-16	IV113S1-16	-	-

2 - WAY VALVE SETS - 10 CM



1,5x2,5mm		1,5x2,7mm		1,5x3mm		3x4,1mm	
UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT
IV210P1	IV210S1	IV211P1	IV211S1	IV212P1	IV212S1	IV213P1	IV213S1
IV220P1	IV220S1	IV221P1	IV221S1	IV222P1	IV222S1	-	-


2 - WAY VALVE SETS - 15 CM



1,5x2,5mm		1,5x2,7mm		1,5x3mm		3x4,1mm	
UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT
IV400P1	IV400S1	IV401P1	IV401S1	IV402P1	IV402S1	-	-
IV410P1-15	IV410S1-15	IV411P1-15	IV411S1-15	IV412P1-15	IV412S1-15	IV413P1-15	IV413S1-15



NEEDLE - FREE CONNECTORS (VALVES) 3 - WAY VALVE SETS - 10 CM



1,5x2,5mm		1,5x2,7mm		1,5x3mm		3x4,1mm	
UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT
IV310P1	IV310S1	IV311P1	IV311S1	IV312P1	IV312S1	IV313P1	IV313S1

3 - WAY VALVE SETS - 15 CM



1,5x2,5mm		1,5x2,7mm		1,5x3mm		3x4,1mm	
UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT	UV PROTECTED	TRANSPARENT
IV510P1	IV510S1	IV511P1	IV511S1	IV512P1	IV512S1	IV513P1	IV513S1



INFUSION ACCESSORIES

As complementary equipment of infusion therapy, there are different types of caps and transfer sets. Safety is given top priority in accessory designs as with the entire product family.



HEPARIN CAP
IV011



TAMPER EVIDENT CAP
SF100



COMBI STOPPER
IV002



BAG SPIKE w. CHECK VALVE & NEEDLE-FREE CONNECTOR
LN112P2



BAG SPIKE w. NEEDLE-FREE CONNECTOR
VA-BSV01



LIQUID TRANSFER SET
ONC509S2



TRANSFER SPIKE
VA-BSV04

OEM


OEM

PRODUCTS



VIAL ADAPTER BODY

Vial adapter bodies with spike, compatible with different brands of needle-free connectors and filters.



REF NO

DESCRIPTION

VA002.20

VIAL ADAPTER BODY WITHOUT FILTER AND VALVE



ERAVERSE

COMPOUNDING DEVICES





ERAVERSE

COMPOUNDING DEVICES

EREVERSE TD – 83

ERAVERSE TC – 84

ERAVERSE MANAGER – 85

TECH SPECS – 86



ERAVERSE TD



Eraverse TD is a drug compounding device that provides safe transfer of drugs and liquids from the source containers to the final container by using closed system transfer sets. Eraverse TD consists of the Eraverse Doser Drug Compounding Module and the Eraverse Touch Control Board. Eraverse Doser Module is available in two types as vertical and horizontal. The module can be designed as single or double channel. Compounding speed can be adjusted according to the type of liquid to be compounded. It provides volumetric compounding.

It is suitable for compounding of TPN, chemotherapy, aseptic IV drug or magistral drugs.



- Provides fast and safe drug compounding.
- Reduces dosing errors.
- Reduces workforce needs.

ERAVERSE TD ERAVERSE TDS

Eraverse TDS is the version which you can verify the weight of the final container with a scale with the help of density information.

- Automatic Filling & Gravimetric Control.
- It helps to eliminate dose errors during preparation process.





ERAVERSE TC



Eraverse TC consists of Eraverse Check Drug Verification Module and Eraverse Touch Control Board. It recognizes the labels of the drugs in seconds and creates a label. It is used before drug preparation process to reduce preparation failures.



- Reduces medication errors.
- Saves time with fast verification.
- Suitable for different dosage forms (Vials, tablets, etc.)

ERAVERSE TC ERAVERSE TCS

Eraverse TCS is the version which you can verify the weight of the final container with a scale with the help of density information.

- Manual Filling & Gravimetric Control.
- It helps to eliminate drug and dosage errors during the preparation process.





ERAVERSE MANAGER SOFTWARE

Eraverse Manager is a drug preparation information management software.

Drug preparation plans can be created by Eraverse Manager and these plans can be forwarded to the compounding modules. Users can be authorized at different security levels. Detailed reports can be prepared by defining drugs, active ingredients or final containers. It is possible to integrate to HIS.



- User-friendly interface.
- Reduces medication errors.
- Automatic documentation.
- Detailed reporting.

Reç Tarifleri

Anahtar sözcükleri buraya

Aktif	Entegrasyon ID	Barkod	Ürün Adı	Ürün Modülü	Kompozit
<input checked="" type="checkbox"/>	1	8699205770041	5-Fluorouracil_360.00mg/1000ml_8bese	5-Fluorouracil	
<input checked="" type="checkbox"/>	2	8699205770034	5-Fluorouracil_1.000.00mg/20.00ml_8bese	5-Fluorouracil	
<input checked="" type="checkbox"/>	3	8699205770005	5-Fluorouracil_3.000.00mg/100.00ml_8bese	5-Fluorouracil	
<input checked="" type="checkbox"/>	4	86992579242	Fluor_5_1.000.00mg/20.00ml_Deve	5-Fluorouracil	
<input checked="" type="checkbox"/>	5	8699402772348	5-Fluor_1.000.00mg/20.00ml_Orko	5-Fluorouracil	
<input checked="" type="checkbox"/>	6	8699402772355	5-Fluor_5.000.00mg/100.00ml_Orko	5-Fluorouracil	
<input checked="" type="checkbox"/>	7	8697942690102	Fluorouracil_250.00mg/5.00ml_Farmako	5-Fluorouracil	
<input checked="" type="checkbox"/>	8	8697942690119	Fluorouracil_500.00mg/10.00ml_Farmako	5-Fluorouracil	
<input checked="" type="checkbox"/>	9	8697942690126	Fluorouracil_1.000.00mg/20.00ml_Farmako	5-Fluorouracil	
<input checked="" type="checkbox"/>	10	8699428770237	Fluorouracil_250.00mg/5.00ml_Kocak	5-Fluorouracil	
<input checked="" type="checkbox"/>	11	8699428770244	Fluorouracil_500.00mg/10.00ml_Kocak	5-Fluorouracil	
<input checked="" type="checkbox"/>	12	8699428770251	Fluorouracil_1.000.00mg/20.00ml_Kocak	5-Fluorouracil	
<input checked="" type="checkbox"/>	13	8699809880530	Zaltrap_200.00mg/8.00ml_Sanofi	Aflibercept	
<input checked="" type="checkbox"/>	14	8699809880537	Zaltrap_330.00mg/4.00ml_Sanofi	Aflibercept	
<input checked="" type="checkbox"/>	15	8699809880545	Zaltrap_500.00mg/10.00ml_Sanofi	Aflibercept	

155

Order Listesi

Anahtar sözcükleri buraya

Seçim	Çinçet	Protokol No	Dozaj No	TC Kimlik No	Hasta Adı Soyadı	Doğum Tarihi	Çalışma Şekli	Order Tarihi
<input checked="" type="checkbox"/>				35433076272	HASAN URŞAR		ONAYLANCIK	25.01.2022
<input checked="" type="checkbox"/>			54321	35433076272	HASAN URŞAR		ONAYLANCIK	09.03.2022
<input checked="" type="checkbox"/>			12345	35433076272	HASAN URŞAR		MANUEL	04.02.2022

Toplam: 3

Uygulandı Uygulama iptal edildi



TECH SPEC

ERAVERSE DOSER - COMPOUNDING MODULE



ERAVERSE DOSER

MODULE TYPE	2 types, horizontal and vertical
NUMBER OF CHANNELS	1 or 2
SYRINGE TYPE	BD 50 ml luer lock syringe
SOURCE CONTAINER TYPE	Vial, bag or syringe
FINAL CONTAINER TYPE	Syringe, bag or elastomeric infusion pump
VOLUME RANGE	0,1 to 9999 mL
INPUT ACCURACY	0.1 mL steps
MAX. COMPOUNDING SPEED	250 mL/min
TRANSFER SET	Different configurations are available based on the source and final container types
DIMENSIONS (WxDxH)	120mm x 250mm x 100mm

- Provides fast and safe drug compounding.
- Reduces dosing errors.
- Reduces workforce needs.

TECH SPEC

ERAVERSE TOUCH - CONTROL BOARD

ERAVERSE TOUCH

DIMENSIONS (WxDxH) 350mm x 200mm x 250mm

SCREEN SIZE 11,6"

APPS

- Eraverse Manager TD
- Eraverse Manager TDS
- Eraverse Manager TC
- Eraverse Manager TCS

- User-friendly interface.
- Warning mechanisms for errors.
- Easy disinfection with disinfectants containing alcohol.





TECH SPEC

ERAVERSE SENSE - WEIGHT VERIFICATION MODULE



ERAVERSE SENSE	
DIMENSIONS (WxDxH)	205mm x 300mm x 100mm
CAPACITY	6000g
READIBILITY	0,1g

- Provides precise and safe dosing.
- Reduces dosing errors.

TECH SPEC

ERAVERSE CHECK - DRUG VERIFICATION MODULE

ERAVERSE CHECK	
DIMENSIONS (WxDxH)	280mm x 190mm x 170mm
DOSAGE FORMS	Vial and medicine packaging box
VALIDATION TIME	Max. 3 sec

- Reduces medication errors.
- Saves time with fast verification.
- Suitable for different dosage forms.
(Vials, tablets, etc.)





COMPOUNDING OF ORAL MIXTURES, IV DRUGS, TPN MIXTURES AND ANTINEOPLASTIC AGENTS BY USING AUTOMATION SYSTEMS REDUCES THE RISKS OF MEDICATION ERRORS.

ONCOERA
Safe for your health

ERASER[®]
MEDICAL

ERASER MEDICAL LLC

KEMALPAŞA MH. PINAR CD. NO:78
PINARBAŞI/BORNOVA, İZMİR - TURKEY

TEL : +90 232 479 80 60

FAX : +90 232 479 80 25

www.eraser.com.tr
info@eraser.com.tr
info@oncoera.com.tr



